

# Real-time Intraoperative Flow Data to Help Surgeons Enhance Outcomes

- Check the Quality of Inflow
- Quantify Flow Improvement
- Confirm Conduit Patency



**transonic**  
THE MEASURE OF BETTER RESULTS.

CE  
2797

# Volume Flow Measurements Add Quantitative Dimension to Vascular Surgery

Vascular surgeons seek to perform the best possible operations tailored to the specific needs of each patient. To this end, intraoperative blood flow measurements provide a quick assessment that may confirm what appears to be an acceptable surgical result, or they may alert the surgeon to potential problems at a time when they can be addressed before closing the patient.

Transonic® volume flow measurements provide on-the-spot feedback about key surgical revascularization/reconstruction components:

- Inflow: checks the quality of flow of a newly constructed fistula or bypass.
- Flow Conduit: quantifies the carrying capacity of the bypass; checks for kinks or twists in a graft.
- Outflow: quantifies flow augmentation and confirms the patency of an anastomosis.

This rapid flow data offers functional information to the vascular surgeon who values real-time data to assess the quality of the surgery. No other flow technology produces flow information so quickly, accurately, and non-intrusively during vascular surgeries as do Transonic® intraoperative Flowmeters.

*"Accurate flow measurements can be of great assistance during vascular reconstructive surgery. The primary aim with these intraoperative measurements is to obtain information on the immediate result of the reconstruction, where a technical failure may jeopardize an otherwise successful operation."*

A Lundell, MD, FACS


*"Transonic Flow-QC provides a measurable improvement in the quality of care you can extend to your patients. With Transonic Flow-QC you can: improve patient outcomes; reduce or delay the need for future interventions and document surgical results."*

T. Wolvos, MD, FACS

*"Intraoperative blood flow measurements obtained at the time of autologous AVF construction can identify fistulas that are unlikely to mature; and therefore, that require immediate revision or abandonment which will ultimately expedite the establishment of a useful access in the HD patient."*

Berman et al, J Vasc Access. 2008;9(4):241-7.

## TRANSIT-TIME ULTRASOUND TECHNOLOGY MEASURES VOLUME FLOW, NOT VELOCITY



Two transducers pass ultrasonic signals, alternately intersecting the vessel in upstream and downstream directions. The difference between the two transit times yields a measure of volume flow.



Transonic Systems Inc. is a global manufacturer of innovative biomedical flow measurement equipment. Founded in 1983, Transonic sells state-of-the-art, transit-time ultrasound devices for surgical, hemodialysis, perfusion, ECMO, and medical device testing applications, and for incorporation into leading edge medical devices.

### USA/Canada

Transonic Systems Inc.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

### Europe

Transonic Europe B.V.  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

### Asia/Pacific

Transonic Asia Inc.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

### Japan

Nipro-Transonic Japan Inc.  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
japan@transonic.com

# Medical Note

## Intraoperative Blood Flow Measurements during Carotid Endarterectomy

Courtesy of Ian Gordon, MD, PhD, Professor, Dept. of Surgery, Univ. of CA, Irvine

### Surgical Approach

No special adjustment in surgical technique is necessary for measurement of blood flow during carotid surgery. The sites on the carotid arteries skeletonized by dissection for vascular clamp placement are identical to those employed for Flowprobe placement. A 10 mm Flowprobe (sometimes an 8 mm Flowprobe is used) is employed for the distal common carotid artery and 6 mm Flowprobes for the origin of the internal and external carotid arteries.

The Flowprobe on the external carotid is placed just distal to the origin of the superior thyroid artery. We perform the pre-endarterectomy flow measurements immediately after administration of systemic heparin. Measurement of all three arteries (Fig. 1) only requires a few minutes.

After the first measurement, vascular clamps are placed and the endarterectomy is performed. After finishing the endarterectomy, removing all clamps, and establishing hemostasis, a repeat flow measurement is performed. Occasionally, a low flow or turbulent flow waveform is detected indicating the presence of a significant stenosis which requires immediate revision. As a consequence, we try not to reverse heparin with protamine until completion of the flow measurements.




Fig 1: Flowprobe Positions: Pre and Post Endarterectomy.

### REFERENCES

Gordon IL, Goldstein LJ, Ceraldi C, Weir P, Tobis J, Stemmer E, "Hemodynamic Effects of Serial Stenoses in an Arterial Model," American College of Surgeons 1991, Surgical Forum, Vol. XLII, pp. 315-319 (Transonic Reference # 192A)

Gordon, I., Weil, J., Williams, R., Wilson, E., "Intraoperative Measurement of Javid Shunt Flow with Transit Time Ultrasound," Annals of Vasc Surg 1994; 8: 571-577. (Transonic Reference # 192A)

Gordon IL, Stemmer EA, Williams RA, Araf M, Wilson SE, "Changes in Internal Carotid Blood Flow after Carotid Endarterectomy Correlates with Preoperative Stenosis," The Am J Surg 1994; 168: 127-130 (Transonic Reference # 364AH).

Gordon IL, Stemmer EA, Wilson SE, "Redistribution of blood flow after carotid endarterectomy," J Vasc Surg 1995; 22: 349-60. (Transonic Reference # 593AH)

TYPICAL FLOWS OBSERVED		
Conduit	Pre-Endarterectomy	Post-Endarterectomy
Common carotid a	287 ± 17	329 ± 18
External carotid a	126 ± 10	104 ± 8
Internal carotid a	135 ± 10	178 ± 11

Carotid Endarterectomy	Probe Size (mm)
Common carotid artery	8, 10
External carotid artery	6
Internal carotid artery	6



6 mm Carotid Flowprobe (-FME)



# Medical Note

## Intraoperative Blood Flow Measurements during Lower Extremity Bypass Grafts

Courtesy of Ian Gordon, MD, PhD, Professor, Dept. of Surgery, Univ. of CA, Irvine

### Introduction

After completing a lower extremity arterial bypass anastomosis, flow is measured immediately before closure of the wound. We do not reverse heparinization.

Three methods are used to measure flow depending on the graft material used. Method A is suitable for saphenous vein or dacron grafts and Methods B and C are useful for PTFE grafts.

### 1: Saphenous Vein/Dacron Bypass Grafts

**Proximal Flow:** To measure flow in a saphenous vein or dacron bypass graft, the Flowprobe is placed just distal to the proximal anastomosis (Fig. 1: A). Flow is documented and a representative flow waveform is generated. An artifact in the flow versus time waveform indicates the presence of a hemodynamically significant stenosis causing turbulence.

**Distal Flow:** The distal anastomosis is similarly assessed for turbulence by placing the Flowprobe on the target vessel for the bypass just distal to the distal anastomosis (Fig. 1: B).

**Bypass Flow:** If there are no technical problems requiring graft revision, we perform our definitive flow measurement by positioning the Flowprobe on the bypass at any convenient position (Fig. 1: C) and measuring flow. We customarily measure flow, first with the graft temporarily clamped to confirm an accurate zero flow. Then graft flow [ $F_{\text{graft}}$ ] is measured with the clamp released, and recorded.

**Resistance:** To calculate resistance of graft flow in the distal run-off vessels, we measure the pressure drop across the graft. A 26 gauge needle, connected to a three-way stopcock, is connected by plastic extension tubing to a sterile pressure transducer (usually the anesthesiologist's radial artery catheter transducer) and is introduced onto the surgical field. The bypass graft is punctured by the needle several centimeters distal to the proximal anastomosis (Fig. 2). Mean pressure is measured first, with the graft open [ $P_{\text{open}}$ ], and then, with a clamp proximal to the needle, occluding the graft [ $P_{\text{clamp}}$ ]. The pressures are recorded. After completing the measurements, the needle is withdrawn and the needle hole is closed with a 6-0 suture.

Flow resistance (R) is calculated as:

$$R = F_{\text{graft}} / (P_{\text{open}} - P_{\text{clamp}})$$




Fig. 1: Flow measurement sites for fem-pop saphenous vein or dacron grafts.

Fig. 2: To calculate resistance, the pressure drop across the bypass graft is measured, first with the clamp on the bypass opened [ $P_{\text{open}}$ ] and then with the clamp closed [ $P_{\text{clamp}}$ ].

Fig. 3: Method B: PTFE bypass grafts: flow is measured proximal to the bypass with the bypass open [ $F_{\text{open}}$ ] and then closed [ $F_{\text{clamp}}$ ]. To measure flow distribution beyond the distal anastomosis, either antegrade or retrograde, flows are measured at points B and E with the bypass open and clamped.

### 2: PTFE Femoropopliteal Bypass

In a typical PTFE femoropopliteal bypass, the graft origin is from the common femoral artery. An 8 or 10 mm Flowprobe is placed on the common femoral artery just proximal to the bypass (Fig. 3, D). Flow in the common femoral artery is then measured with the graft open [ $F_{\text{open}}$ ] or clamped [ $F_{\text{clamp}}$ ]. Net graft flow is equal to ( $F_{\text{open}} - F_{\text{clamp}}$ ). Resistance is measured as above.

Method B is also used whenever placement of the Flowprobe on both sides of the distal anastomosis is difficult due to poor exposure of the distal vessel receiving the bypass.

**Note:** Dacron grafts, rarely used distal to the femoral artery, allow direct measurement of flow by transit-time Flowprobes. Expanded PTFE grafts cannot be studied immediately by a Flowprobe placed on the graft, as air trapped in the graft interstices interferes with ultrasound transmission, and an accurate measurement is not possible until this gas is expelled.

# Intraoperative Blood Flow Measurements during Lower Extremity Bypass Grafts Cont.

## 3: Antegrade & Retrograde Flow Distribution

Method C is used to measure the distribution of flow beyond the distal anastomosis in retrograde and antegrade directions. A 4 or 6 mm probe is placed on the target vessel just distal to the distal anastomosis [position B].

**Antegrade Flow:** Antegrade flow is measured with the bypass graft open [ $F_{ao}$ ] and clamped [ $F_{ac}$ ].

**Retrograde Flow:** To measure retrograde flow, the probe is placed on the target vessel [position E] just proximal to the distal anastomosis and flow is measured with the graft open [ $F_{ro}$ ] and clamped [ $F_{rc}$ ]. Observe the direction of blood flow carefully observed and negative and positive signs correctly employed to accurately measure flow. Net graft flow is calculated as  $(F_{ao} - F_{ac}) + (F_{ro} - F_{rc})$ . Resistance is measured as in Method A.

Note: Pulse is a manifestation of pressure, not flow, so an occluded graft may still have a distinct pulse.




Fig. 4: These waveforms demonstrate a case in which a vein graft initially had a stenosis near its origin. The upper waveform shows flow (45 mL/min) before correcting the stenosis. The lower waveform shows flow (90 mL/min) along with a representative flow pattern after correcting the stenosis.

## Equipment Needs

BYPASS	SIZE	PROBE SERIES
Profunda femoris	8	-FMV, -AU
Common femoral	8, 10	-FMV, -AU
Popliteal	4, 6	-FMV, -AU
Tibial	3, 4	-FMV, -AU



2 mm - 10 mm Vascular Handle Flowprobe (-FMV-Series)



Transonic Systems Inc. is a global manufacturer of innovative biomedical flow measurement equipment. Founded in 1983, Transonic sells state-of-the-art, transit-time ultrasound devices for surgical, hemodialysis, perfusion, ECMO, and medical device testing applications, and for incorporation into leading edge medical devices..

### AMERICAS

Transonic Systems Inc.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

### EUROPE

Transonic Europe BV  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

### ASIA/PACIFIC

Transonic Asia Inc.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

### JAPAN

Nipro-Transonic Japan Inc.  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
japan@transonic.com

# Vascular Surgery References

## AV Access

### Creation

Johnson CP, Zhu Y, Matt C, Pelz C, Roza AM, Adams MB, "Prognostic Value of Intraoperative Blood Flow Measurements in Vascular Access Surgery," *Surgery* 1998; 124: 729-38. (1504AH)

Berman SS, Mendoza B, Westerband A, Quick RC, "Predicting Arteriovenous Fistula Maturation with Intraoperative Blood Flow Measurements," *Journal of Vascular Access* 2008; 9(4): 241-7. (7710AH)

Won T, Jang JW, Lee S, Han JJ, Park YS, Ahn JH, "Effects of intraoperative blood flow on the early patency of radiocephalic fistulas," *Ann Vasc Surg* 2000; 14(5): 468-72. (2411AH)

Lin CH, Chua CH, Chiang SS, Liou JY, Hung HF, Chang CH., "Correlation of Intraoperative Blood Flow Measurement with Autogenous Arteriovenous Fistula Outcome." *J Vasc Surg.* 2008; 48: 167-72 (7637AH)

### Flow Reduction

van Hoek F, Scheltinga, MR, Kouwenberg, I, Moret, KEM. Beerenhout, CH, Tordoir, JHM. "Steal in hemodialysis patients depends on type of vascular access," *Eur J Vasc Endovasc Surg* 2006; 32: 710-717. (HD7573A)

van Hoek F, Scheltinga MR, Luirink M, Raaijmakers LCJ, van Pul C, Beerenhout CH, "Access flow, venous saturation and digital pressures in hemodialysis," *J Vasc Surg* 2007;45: 968-73. (7572AH)

van Hoek F, Scheltinga MR, Luirink M, Pasmans H, Beerenhout CH, "Banding of hemodialysis access causing hand ischemia or cardiac overload," *Seminars in Dialysis* 2009; 22: 204-208. (9719AH)

Scheltinga MR, van Hoek F, Bruyninckx CMA, "Surgical Banding for Refractory Hemodialysis Access-Induced Distal Ischemia (HAIDI)," *J Vasc Acc* 2009;10: 43-49. (9720AH)

van Hoek F, In: *Vascular Access*. J. Tordoir (ed). "Banding for high flow hemodialysis access (HFA)," *Minerva*, Turin, 2009, 141-150.

Chemla ES, Tang VC, Eymann SA, "Intraoperative Flow Measurements Are Helpful in the Treatment of High-Inflow Steal Syndrome on a Predialysis Patient with a Brachiocephalic Fistula: A Case Report," *Ann Vasc Surg.* 2007. (7402AH)

Zanow J, Petzold K, Petzold M, Kruege U, Scholz H, Flow reduction in high-flow arteriovenous access using intraoperative flow monitoring, *J Vasc Surg* 2006 Dec;44(6):1273-8. (7410AH)

### Other

Jendrisak MD, Anderson CB, "Vascular Access in Patients with Arterial Insufficiency," *Ann Surg* 1990; 212(2): 187-193. (185AH)

Kelber J, Deimez JA, Windus DW, "Factors Affecting Delivery of High-Efficiency Dialysis Using Temporary Vascular Access," *Am J Kid Dis* 1993; 22(1): 24-29. (365AH)

Wong V, Taylor J, Selvakumar S, How TV, Bakran A., "Factors associated with early failure of arteriovenous fistulae for haemodialysis access." *Eur J Vasc Endovasc Surg* 1996; 12(2): 207-13. (6943AH)

Ikizler TA, Pupim LB, Brouillette JR, Levenhagen DK, Farmer K, Hakim RM, Flakoll PJ, "Hemodialysis stimulates muscle and whole body protein loss and alters substrate oxidation," *AJP* 2002; 282(1): E107-16. (2636AH)

## Peripheral Vascular

Bosma J, Minnee RC, Erdogan D, Wisselink W, Vahl AC, "Transit-time Volume Flow Measurements in Autogenous Femorodistal Bypass Surgery for Intraoperative Quality Control," *Vascular* 20120; 18(6):344-9. (8091AH)

Desai TR, Meyerson SL, Skelly CL, MacKenzie KS, Bassiouny HS, Katz D, McKinsey JF, Gewertz BL, Schwartz LB, "Patency and limb salvage after infrainguinal bypass with severely compromised ("blind") outflow," *Arch Surg* 2001; 136(6): 635-42. (7551AHR)

Curi MA, Skelly CL, Quint C, Meyerson SL, Farmer AJ, Shakur UM, Loth F, Schwartz LB. "Longitudinal impedance is independent of outflow resistance," *J Surg Res* 2002; 108(2): 191-7. (7240AH)

Lorenzetti F, Tukianen E, Alback A, Kallio M, Asko-Seijavaara S, Lepantalo M, "Blood Flow in a Pedal Bypass Combined with a Free Muscle Flap," *Eur J Endovasc Surg* 2001; 22: 161-164. (2709AHM)

Meyerson SL, Skelly CL, Curi MA, Desai TR, Katz D, Bassiouny HS, McKinsey JF, Gewertz BL, Schwartz LB, "Long-term results justify autogenous infrainguinal bypass grafting in patients with end-stage renal failure," *J Vasc Surg* 2001; 34(1): 27-33. (7293AH)

Stirnemann P, Ris HB, Do D, Hammerli R, "Intraoperative flow measurement of distal runoff: a valid predictor of outcome of infrainguinal bypass surgery," *Eur J Surg.* 1994; 160(8): 431-6. (7249AH)

Tukiainen E, Laurila K, Kallio M, Lorenzetti F, Kantonen I, Lepantalo M, "Internal arteriovenous fistula within a radial forearm flap—a novel technique to increase femorodistal bypass graft flow to the diabetic foot and flap covering ischaemic tissue loss," *Eur J Vasc Endovasc Surg* 2006; 31(4): 423-30. (7251AHM)

Schwierz T, Harnoncourt F, Havlicek W, Tomaselli F, Fuggerm R, "Interpretation of the results of Doppler ultrasound flow volume measurements of infrainguinal vein bypasses," *Eur J Vasc Endovasc Surg*, 2005; 29(5): 452-6. (7246AHM)

Schwierz T, Pricop T, Ebner C, Havlicek W, Schmoller F, Nesser HJ, Bohmig HJ, Fugger R. "The evaluation of run-off prior to infra-inguinal bypass reconstruction - a modified scoring system based on flow measurement," *Eur J Vasc Endovasc Surg* 2003; 26(1): 52-8. (7245AHM)

Fisher RK, How TV, Bakran A, Brennan JA, Harris PL, "Outflow distribution at the distal anastomosis of infra-inguinal bypass grafts," *J Biomech* 2004; 37(3): 417-20. (2943AH)

Meyerson SL, Moawad J, Loth, F, Skelly, CL, Bassiouny, HS, McKinsey, JF, Sewertz, BL, Schwartz, LB, "Effective Hemodynamic Diameter: An Intrinsic Property of Vein Grafts with Predictive Value for Patency," *J Vasc Surg* 2000; 31: 910-917. (1503AH)

Schwartz, LB, Belkin, M, Donaldson, MC, Knox, JB, Craig, DM, Moawad J, McKinsey JF, Piano G, Bassiouny HS, Whittemore AD, "Validation of a New and Specific Intraoperative Measurement of Vein Graft Resistance" *J Vasc Surg* 1997; 25: 1033-1043. (925AH)

Lundell A, Begqvist D, "Intraoperative Flow Measurements in Vascular Reconstruction," *Annales Chirurgiae Gynaecologiae* 1992; 81(2): 187-191. (3G)

# Vascular Surgery References Cont.

## Peripheral Vascular cont.

Lundell, A, Nyborg, K, "Do Residual Arteriovenous Fistulae after in Situ Saphenous Vein Bypass Grafting Influence Patency?" J Vasc Surg 1999; 30: 99-105. (1903AH)

Lundell, A, Bergqvist, D, "Prediction of Early Graft Occlusion in Femoropopliteal and Femorodistal Reconstruction by Measurement of Volume Flow with a Transit Time Flowmeter and Calculation of Peripheral Resistance," Eur J Vasc Surg 1993; 7: 704-9. (293AH)

Pedersen G, Laxdal E, Amundsen SR, Dregelid E, Jonung T, Nyheim T, Aune S, Flow measurement before and after papaverine injection in above-knee prosthetic femoro-popliteal bypass," J Vasc Surg 2006; 43(4): 729-34. (7136AHM)

Alback A, Roth WD, Ihlberg L, Biancari F, Kepantalo M, "Preoperative angiographic score and intraoperative flow as predictors of the mid-term patency of infrapopliteal bypass grafts," Eur J Vasc Endovasc Surg 2000; 20(5): 447-53. (2062AHM)

Ihlberg LHM, Alback NA, Lassila R., Lepantalo M, "Intraoperative Flow Predicts the Development of Stenosis in Infringuinal Vein Grafts," J Vasc Surg 2001; 34: 269-76. (2064AHM)

Blankensteijn JD, Gertler JP, Brewster DC, Cambria RP, LaMurgia GM, Abbott WM, "Intraoperative Determinants of Infringuinal Bypass Graft Patency: A Prospective Study," Eur J Endovasc Surg 1995; 9: 375-382. (518AH)

Belkin M, Schwartz LB, Donaldson MC, Mannick, JA Whittemore, AD, "Hemodynamic Impact of Vein Graft Stenoses and Their Prediction in the Vascular Laboratory," J Vasc Surg 1997; 25: 1016-1022. (924AH)

Verhelst R, Bruneau M, Nicolas A-L, Frangi R, Khoury GE, Noirhomme P, Dion R, "Popliteal-to-Distal Bypass Grafts for Limb Salvage," Ann Vasc Surg 1997; 11: 505-509. (1043AH)

## Carotid

Hosoda K, Kawaguchi T, Shibata Y, Kamei M, Kidoguchi K, Koyama J, Fujita S, Tamaki N, "Cerebral vaso-reactivity and internal carotid artery flow help to identify patients at risk for hyperperfusion after carotid endarterectomy," Stroke 2001; 32(7): 1567-73. (6995AH)

Eckstein HH, Eichbaum M, Klemm K, Doerfler A, Ringleb P, Bruckner T, Allenberg JR, "Improvement of Carotid Blood Flow after Carotid Endarterectomy—Evaluation using Intraoperative Ultrasound Flow Measurement," Eur J Vasc Endovasc Surg 2003; 25: 168-174. (2579AH)

Gordon IL, Stemmer EA, Wilson SE, "Redistribution of Blood Flow after Carotid Endarterectomy," J Vasc Surg 1995; 22: 349-360. (593AH)

Gordon IL, Weil J, Williams RA, Wilson SE, "Intraoperative Measurement of Javid Shunt Flow with Transit-Time Ultrasound," Ann Vasc Surg 1994; 8(6): 571-577. (371AH)

Gordon IL, Stemmer EA, Williams RA, Arafat M, Wilson SE, "Changes in Internal Carotid Blood Flow after Carotid Endarterectomy Correlate with Preoperative Stenosis," The Am J Surg 1994; 168: 27-130. (364AH)

Gordon IL, Goldstein LJ, Ceraldi C, Weir P, Tobis J, Stemmer E, "Hemodynamic Effects of Serial Stenoses in An Arterial Model," Am College Surgeons, Surgical Forum 1991; XLII: 315-319. (192AH)

Aleksic M, Heckenkamp J, Gawenda M, Brunkwall J, "Pulsatility index determination by flowmeter measurement: a new indicator for vascular resistance?" Eur Surg Res 2004; 36(6): 345-9. (3018AHM)

McCullough JN, Zhang N, Reich DL, Juvonen TS, Klein JJ, Spielvogel D, Ergin MA, Griep RB, "Cerebral metabolic suppression during hypothermic circulatory arrest in humans," Ann Thorac Surg 1999; 67(6): 1895-9; discussion 1919-21. (7191AH)

Tulleken CA, Verdaasdonk RM, Beck HJM, "Nonocclusive Excimer Laser-Assisted End-to-Side Anastomosis," Ann Thorac Surg 1997; 63: S138-142. (1037AH)

Tulleken CA, van der Zan A, Van Rooij WJ, Ramos LM, "High-flow Bypass Using Nonocclusive Excimer Laser-assisted End-to-side Anastomosis of the External Carotid Artery to the P1 Segment of the Posterior Cerebral Artery Via the Sylvian Route. Technical Note," J Neurosurg 1998; 88(5): 925-927. (1025AH)

Tulleken CA, van der Zwan A, Kappelle LJ, "High-flow transcranial bypass for prevention of brain ischemia," Ned Tijdschr Geneesk 1999; 143(45): 2281-5. (3028AH)

## Other

Murakawa K, Ishimoto E, Noma K, Ishida K, Nishijima M, Izumi R, "Circulatory Effects of Stellate Ganglion Block in Idiopathic Facial Palsy," Masui 1994; 43(3): 356-60, 1994. (1005AH)

Schmidli J, Heller G, Englberger L, Eckstein FS, Carrel TP, "First clinical experience with an anastomotic device to facilitate aortomesenteric saphenous vein bypass," J Vasc Surg 2002; 36(4): 859-62. (2760AH)

Herz I, Mohr R, Aviram G, Loberman D, Locker C, Ben-Gal Y, Uretzky G, "The right internal thoracic artery and right gastroepiploic artery: alternative sites for proximal anastomosis in patients with atherosclerotic calcified aorta," Heart Surg Forum - 2004; 7(5): E481-4. (6942AH)



Transonic Systems Inc. is a global manufacturer of innovative biomedical flow measurement equipment. Founded in 1983, Transonic sells state-of-the-art, transit-time ultrasound devices for surgical, hemodialysis, perfusion, ECMO, and medical device testing applications, and for incorporation into leading edge medical devices.

### AMERICAS

Transonic Systems Inc.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

### EUROPE

Transonic Europe B.V.  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

### ASIA/PACIFIC

Transonic Asia Inc.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

### JAPAN

Transonic Japan Inc.  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
japan@transonic.com