

Flow-assisted Surgical Techniques and Notes*

Arterial EC-IC Bypass Surgery Protocol

Drawn from the clinical expertise of FT Charbel, MD, FACS, S Amin-Hanjani, Univ. of IL at Chicago, Chicago, IL *et al.*¹⁻²⁴

*Flow-Assisted Surgical Techniques (“F•A•S•T”) and Protocols are drawn from surgical experiences by transit-time flow measurement users and passed along by Transonic for educational purposes. They are not intended to be used as sole basis for diagnosis. Clinical interpretation of each patient’s individual case is required.

Introduction³

An extracranial to intracranial (EC-IC) bypass is used during cerebrovascular surgery:

- 1) to augment flow for occlusive cerebrovascular disease (i.e., Moyamoya)
- 2) to replace flow during aneurysm clipping surgery when an aneurysm is trapped and a parent vessel (i.e., internal carotid) has to be occluded and sacrificed.

Flow Augmentation for Occlusive Cerebrovascular Disease¹⁻⁵

In 2005, Drs. Fady Charbel and Sepideh Amin-Hanjani introduced the concept of a Cut Flow Index to evaluate the quality of an EC-IC bypass used to enhance flow during cerebral ischemia. Briefly, the free flow of the donor extracranial artery intended for use as a bypass is measured. Once the bypass is constructed, the bypass flow of the donor artery is measured. The ratio of bypass flow to free flow is the Cut Flow Index. A value greater than 0.5 indicates that the bypass should be viable.

Flow Replacement during Aneurysm Clipping Surgery^{1,6}

Dr. Amin-Hanjani developed a strategy to assess the adequacy of an STA-or an occipital artery bypass to replace flow when an aneurysm has to be trapped and a parent vessel sacrificed.

Flow Deficit Determined

Flow in the artery or territory distal to the aneurysm is measured and recorded. The vessel to be sacrificed is temporarily occluded and flow is again measured in the distal artery or territory. The difference between the two flows represents the amount of flow deficit that can be expected if the parent vessel is sacrificed. This is the flow that the bypass will have to replace.

Free Flow Determined

The “Free” or “Cut Flow” of the intended bypass is then measured. This Cut Flow value is compared to Deficit Flow. If the Cut Flow value equals or exceeds the potential flow deficit, the EC-IC bypass is completed and the vessel can be sacrificed with reasonable assurance that the bypass flow will compensate for the flow deficit from the sacrificed parent vessel.

Example: STA to M3 Bypass (ICA Aneurysm Clipped, Trapped and ICA Sacrificed)^{1,6}

- | | |
|--|------------------|
| 1) M1 baseline flow measured | 70mL/min |
| 2) M1 flow measured with ICA temporarily occluded | <u>50 mL/min</u> |
| 3) Anticipated Flow Deficit Calculated
(if aneurysm trapped and parent vessel sacrificed) | 20 mL/min |
| 4) STA Cut Flow measured
(STA bypass should be able to supply the flow deficit) | 44 mL/min |
| 5) STA Bypass to M3 completed; aneurysm clipped and trapped | |
| 6) STA Bypass Graft Flow measured
(bypass flow can compensate for anticipated flow deficit) | 24 mL/min |

Dr. Amin-Hanjani reported that this selective strategy allows the surgeon to:

- 1) Assess the adequacy of a bypass before completing its construction
- 2) Select the best match for a bypass
- 3) Evaluate the bypass immediately

Flow-assisted Surgical Techniques and Notes


Arterial EC-IC Bypass Surgery Protocol cont.

When a surgeon selects an arterial extracranial-intracranial (EC-IC) Bypass to preserve flow during aneurysm clipping or trapping surgery, Charbel Probes® assess the adequacy of flow(s) during and after construction of the bypass.¹⁻⁵

Extracranial Donor Artery¹⁻⁷

1. Choose an appropriate-sized Charbel Probe® for the donor (STA) artery.

PROBE SIZE	VESSEL RANGE, OUTER DIAMETER
1.5 mm	1.1 - 1.6 mm
2 mm	1.6 - 2.4 mm
3 mm	2.6 - 3.8 mm



2. Measure baseline flow in the donor artery. Record on the EC-IC Bypass Record.
3. Cut the donor extracranial artery and measure the artery's "Free Flow" by allowing the cut distal end to bleed freely for 15-20 seconds (Fig. 1). This free flow or "Cut Flow" is the amount of flow at zero resistance or the "carrying" capacity or maximum flow the artery can deliver. Record flow on the EC-IC Bypass Record.⁵





Fig. 1: Measurement of "cut" flow in donor artery.

Intracranial Recipient Artery¹⁻⁷

1. Choose an appropriate size Charbel Probe® for recipient artery.

PROBE SIZE	VESSEL RANGE, OUTER DIAMETER
1.5 mm	1.1 - 1.6 mm
2 mm	1.6 - 2.4 mm
3 mm	2.6 - 3.8 mm



2. Measure and record baseline flow in recipient intracranial artery distal to target anastomotic site.
3. Re-measure, record baseline flow in recipient intracranial artery distal to target anastomotic site with vessel to be sacrificed occluded.

4. Calculate anticipated flow deficit by subtracting flow with vessel occluded from baseline recipient arterial flow.




Fig. 2: Flow measurement of bypass after anastomosis to recipient artery.

Construct EC-IC Bypass¹⁻⁷

1. Anastomose the extracranial bypass to the recipient arterial vessel.
2. Measure post-bypass flow in the donor artery (Fig. 2). Record flow.
3. Calculate the Cut Flow Index (CFI) by dividing the Post-Bypass Flow by the Free or Cut Flow (Fig. 1).
 - If post-bypass flow exceeds 50% of (CFI > 0.5), the bypass can be considered successful.
 - If bypass flow is below 50% of free flow (CFI < 0.5), examine bypass for kinks, analyze recipient bed.

Cerebrovascular EC-IC Bypass References

1. Amin-Hanjani S, Alaraj A, Charbel FT, "Flow replacement bypass for aneurysms: decision-making using intraoperative blood flow measurements. *Acta Neurochir (Wien)*. 2010 Jun; 152(6):1021-32; discussion 1032. (Transonic Reference # 7940AH)
2. Amin-Hanjani S, "Cerebral revascularization: extracranial-intracranial bypass," *J Neurosurg Sci*. 2011 Jun; 55(2): 107-16. (Transonic Reference # 10097AH) 6
3. Cerebrovascular Surgery Handbook NS-59-hb, Rev F, 2018
4. Measuring PeriFlowprobe(CV-180-mn) RevA2018US1tr
5. Amin-Hanjani S, Charbel FT *et al*, "The Cut Flow Index: An Intraoperative Predictor of the Success of EC-IC Bypass for Occlusive Cerebrovascular Disease," *Neurosurgery* 2005; 56: 75-85. (Transonic Reference # 2922AH)
6. Amin-Hanjani S, "EC-IC Bypass for Aneurysm: Decision Making Using Intraoperative Flow Measurements," 8th International Conference on Cerebrovascular Surgery (Tapei, 2006).
7. Charbel FT, Meglio G, Amin-Hanjani S, "Superficial Temporal Artery - Middle Cerebral Artery Bypass," *Neurosurgery* 2005; 56: 186-90. (Transonic Reference # 2921AH)
8. Nossek E, Langer DJ, "How I do it: combined direct (STA-MCA) and indirect (EDAS) EC-IC bypass," *Acta Neurochir (Wien)*. 2014; 156(11): 2079-84. (Transonic Ref # 10076AH)
9. Kaku Y, Funatsu N, Tsujimoto M, Yamashita K, Kokuzawa J, "STA-MCA/STA-PCA Bypass Using Short Interposition Vein Graft," *Acta Neurochir Suppl*. 2014; 119: 79-82. (Transonic Reference # 10085AH)
10. van der Zwan A, "How I Do It:" Non-occlusive High Flow Bypass Surgery, *Acta Neurochir Suppl*. 2014; 119: 71-6. (Transonic Reference # 10086AH)
11. Amin-Hanjani S *et al*, "Resolution of bilateral moyamoya associated collateral vessel aneurysms: Rationale for endovascular versus surgical intervention," *Surg Neurol Int*. 2014 Jun 19; 5(Suppl 4): S155-60. (Transonic Reference # 10095AH)
12. Amin-Hanjani S, Charbel FT *et al*, "Combined Direct and Indirect Bypass for Moyamoya: Quantitative Assessment of Direct Bypass Flow Over Time," *Neurosurgery* 2013; 73(6): 962-8. (Transonic Reference # 9835AH)
13. Lee M, Steinberg GK *et al*, "Intraoperative blood flow analysis of direct revascularization procedures in patients with moyamoya disease. *J Cereb Blood Flow Metab*. 2011 Jan; 31(1): 262-74. (Transonic Reference # 7969AH)
14. Kim JY *et al*, "Changes in Bypass Flow during Temporary Occlusion of Unused Branch of Superficial Temporal Artery," *J Korean Neurosurg Soc*. 2010; 48: 105-108. (Transonic Reference # 8006AH)
15. Amin-Hanjani S, "Diagnosis and neurosurgical treatment of intracranial vascular occlusive syndromes." *Curr Treat Options Cardiovasc Med* 2009 Jun; 11(3): 212-20. (Transonic Reference # 7821AH)
16. Guzman R, Steinberg GK, "Direct bypass techniques for the treatment of pediatric moyamoya disease," *Neurosurg Clin N Am*. 2010 Jul; 21(3): 565-73. (Transonic Reference # 8010AH)
17. Lee M, Steinberg GK *et al*, Quantitative hemodynamic studies in moyamoya disease: a review, *Neurosurg Focus*. 2009 Apr; 26(4):E5. (Transonic Reference # 7762AH)
18. Amin-Hanjani S, Charbel FT, "Is extracranial-intracranial bypass surgery effective in certain patients?" *Neurosurg Clin N Am*. 2008; 19(3): 477-87, vi-vii (Transonic Reference # 7703AH)
19. Ashley WW, Charbel FT *et al*, "Flow-assisted surgical cerebral revascularization," *Neurosurg Focus*. 2008; 24(2): E20. (Transonic Reference # 7587AH)
20. Alaraj A *et al*, "The superficial temporal artery trunk as a donor vessel in cerebral revascularization: benefits and pitfalls," *Neurosurg Focus*. 2008; 24(2): E7. (Transonic Reference # 7588AH)
21. Amin-Hanjani S, Charbel FT, "Flow-assisted surgical technique in cerebrovascular surgery," *Surg Neurol*. 2007; 68 Suppl 1:S4-11. (Transonic Reference # 7560AHR)
22. Nossek E, Langer DJ *et al*, "Infratemporal approach for subcranial-intracranial (SC-IC) bypass," *Neurosurg* 2014; 75(1): 87-95 (Transonic Ref # 10061AH)
23. Nossek E, Chalif DJ, Dehdashti AR, "How I do it: occipital artery to posterior inferior cerebellar artery bypass," *Acta Neurochir (Wien)*. 2014 May; 156(5): 971-5. (Transonic Ref # 10077AH).
24. Dunn GP, Gerrard JL, Jho DH, Ogilvy CS, "Surgical treatment of a large fusiform distal anterior cerebral artery aneurysm with *In Situ* end-to-side A3-A3 bypass graft and aneurysm trapping: case report and review of the literature," *Neurosurgery* 2011; 68(2): E587-91. (Transonic Reference # 10084AH)
25. AU-QRG-Optima-EN Rev E