

Transplant Medical Note

Portal Vein Blood Flow Measurement during Auto Islet Cell Transplantation after Pancreatectomy

Flow Measurement during Islet Infusion

Excising a diseased pancreas removes not only pancreatic cells that produce digestive enzymes but also islet of Langerhans cells that produce insulin to control blood sugar. Without insulin a patient becomes diabetic and requires lifelong use of insulin to control blood sugars.

Auto islet cell transplantation takes these islet of Langerhans cells from the pancreas and transplants them to the liver to reduce the diabetic risk. To do this, the removed pancreas is processed to isolate the insulin-producing islets of Langerhans cells. The isolated cells are suspended in a solution and are then slowly infused through the splenic vein back into the patient's liver where it is hoped that they will implant, grow and produce insulin to metabolize sugar.

Typically, 800 - 1500 cc of solution is infused into the portal vein distal to the splenic vein (Fig. 2) over an extended period of time. The team may elect to infuse a small amount over 5 minutes and allow the patient to recover before resuming the infusion. Blood pressure and flow are monitored continuously and for ten minutes after the infusion is completed (Fig. 1).

Flow Measurement during Islet Infusion

Surgeons measure portal venous flow during islet cell infusion to detect any sudden decrease in flow that may foreshadow a problem with the infusion. A 10 mm to 14 mm Perivascular Flowprobe is placed on the portal vein and flow is measured continuously. The Flowprobe is chosen to comfortably encompass - but not constrict - the portal vein. If needed, saline can be used to provide acoustic contact between the vein and Flowprobe. Readings stabilize within 1-2 minutes. Wide fluctuation of measurements may indicate improper positioning of the Flowprobe with poor alignment or fat within the ultrasonic sensing window. Repositioning can normally correct this problem.

Discussion

In this high stakes auto islet cell transplantation procedure, Flowprobes provide a continuous volumetric measure of portal vein flow to inform the surgeon about the safety, fluidity and success of auto islet cell transplantation.




Fig. 1: Steps: Auto Islet Cell Transplantation
 a) Removal of pancreas (pancreatectomy)
 b) Isolation of Islet cells from removed pancreas
 c) Islet cells placed in Infusion bag with solution
 d) Islet cells infused into splenic vein
 e) Islet cells implant in liver




Fig. 2: Enlarged view of islet cell infusion into the splenic/portal venous system.

Portal Vein Blood Flow Measurement during Auto Islet Cell Transplantation after Pancreatectomy Cont.

Equipment Needs



HT354 Single-channel Optima Flowmeter. Acquire precise actual flow measurement quickly, easily and cost effectively.



8 mm to 14 mm FMV Vascular Handle Flowprobes are recommended for portal venous flow measurements during islet cell infusion.



Confidence Flowprobe®

Confidence Flowprobes® provide highly accurate measurements in vessels with fluctuating flows such as the portal vein. The Probes may be left in place for extended measurements and then easily removed via a ring attached to the pliable liner that cushions and protects the vessel.

References

Sutherland DE, Radosevich DM, Bellin MD, Hering BJ, Beilman GJ, Dunn TB, Chinnakotla S, Vickers SM, Bland B, Balamurugan AN, Freeman ML, Pruett TL. "Total pancreatectomy and islet autotransplantation for chronic pancreatitis," J Am Coll Surg. 2012 Apr;214(4):409-24.

Bramis K, Gordon-Weeks AN, Friend PJ, Bastin E, Burls A, Silva MA, Dennison AR, "Systematic review of total pancreatectomy and islet autotransplantation for chronic pancreatitis." Br J Surg. 2012 Jun;99(6):761-6.

http://www.hopkinsmedicine.org/transplant/programs/auto_islet/description.html#total_pancreatectomy

Henderson JM et al, "Hemodynamics During Liver Transplantation: The Interactions Between Cardiac Output and Portal Venous and Hepatic Arterial Flows," Hepatology 1992; 16(3): 715-718.

Henderson JM et al, "Volumetric and Functional Liver Blood Flow Are Both Increased in the Human Transplanted Liver," J Hepatology 1993; 17: 204-207.

Troisi R, de Hemptinne B, "Clinical Relevance of Adapting Portal Vein Flow in Living Donor Liver Transplantation in Adult Patients, Liver Transplantation 2004; 9(9): S36-S41.



www.transonic.com

Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

AMERICAS

Transonic Systems Inc.
34 Dutch Mill Rd
Ithaca, NY 14850
U.S.A.
Tel: +1 607-257-5300
Fax: +1 607-257-7256
support@transonic.com

EUROPE

Transonic Europe B.V.
Business Park Stein 205
6181 MB Elsloo
The Netherlands
Tel: +31 43-407-7200
Fax: +31 43-407-7201
europe@transonic.com

ASIA/PACIFIC

Transonic Asia Inc.
6F-3 No 5 Hangsiang Rd
Dayuan, Taoyuan County
33747 Taiwan, R.O.C.
Tel: +886 3399-5806
Fax: +886 3399-5805
support@transonicasia.com

JAPAN

Transonic Japan Inc.
KS Bldg 201, 735-4 Kita-Akitsu
Tokorozawa Saitama
359-0038 Japan
Tel: +81 04-2946-8541
Fax: +81 04-2946-8542
info@transonic.jp