

# Medical Note

## Intraoperative Blood Flow Measurement during Adult Orthotopic Liver Transplantation

Courtesy of J. Michael Henderson, M.D., F.A.C.S.

### Introduction

Abnormal hepatic hemodynamics and physiology in the transplanted liver pose continuing challenges for the surgeon. A practical method for measuring two of these hemodynamic parameters, portal venous and hepatic arterial flows, is by intraoperative flow measurements. Transit-time ultrasound technology is well suited to measure these flows. Flowprobes are easily applied and do not have to be applied tightly to vessels; they simply encompass the vessel.

### Surgical Approach

Measurement of portal venous and hepatic arterial flows can be easily done at the completion of orthotopic liver transplantation using Transonic Flowprobes. Following completion of the vascular anastomoses, the new liver is reperfused, and hemostasis achieved. Prior to biliary reconstruction, the Flowprobes are placed on the reconstructed portal vein and hepatic artery.

The Probes are chosen to comfortably encompass - but not constrict - the vessels, and are placed such that extraneous tissue is excluded. The field is then immersed in saline which serves as a good acoustic contact with the vessels. Readings stabilize rapidly, usually within 1-2 minutes, and in stable patients fluctuate less than  $\pm 10\%$  when left in situ for 10-15 minutes. If there is wider fluctuation, this usually indicates improper positioning of the Flowprobes with poor alignment or extraneous tissue, and can normally be corrected by repositioning. Arterial flow readings are meaningful over a brief snapshot period. Venous flow exhibits a far slower rhythm, dictated by events such as gastric motility. A one-to-five minute observation period is often adequate.

### Discussion

Combined portal venous and hepatic artery flow are usually 15 - 25% of cardiac output. Of clinical importance is hepatic artery patency and flow, as survival of the graft depends on this. Flowprobes provide a volumetric measure of hepatic artery flow, and when this is low can be used to determine if there is a fixed anatomic limitation to flow or a physiologic limitation. For example, in a patient with a cardiac output of 10 L/min, portal flow of 2000 ml/min and hepatic artery flow of 75 ml/min, reduction of portal flow to 1000 ml/min resulted in a hepatic artery flow increase to 125 ml/min. Thus, the low basal hepatic artery flow resulted from a high physiologic resistance rather than a fixed, potentially surgically correctable low inflow. This kind of data can be collected on the flowmeter's strip chart recorder for a permanent record.

The information obtained with these transit-time ultrasound Flowprobes is often at variance with "clinical impression." A transplant with obstructed hepatic artery may show a strong pressure pulse on the artery, and a healthy organ color due to its venous perfusion. Accurate information on volumetric flow at the time of operation can either be reassuring, or may indicate an unexpected problem which can be fixed at this time.

# Intraoperative Blood Flow Measurements

## Flow - Assisted Liver Transplantation

LIVER HEMODYNAMICS		
TRANSPLANTED LIVER (N = 34) <sup>1</sup>		
Vessel	Flow: Mean ± SD (L/mm)	Range
Total Liver	2.091 ± .932	.570 - 4.540
Portal vein	1.808 ± .929	.300 - 4.500
Hepatic artery	0.268 ± .157	0.30 - 0.675

In a procedure such as liver transplant, where the stakes are high, this technology can be a useful adjunct in operative decision. Subsequent studies have identified the following intraoperative flow indices related to poor outcomes:

- Poor outcome is associated with graft hyperfusion. Recipient portal venous flow in the recipient should be lowered when graft to recipient body weight ratio (GRBWR) < 0.8 is accompanied by portal inflow of > 250 mL/min/100g graft weight.<sup>3</sup>
- Hepatic arterial flow < 100 mL/min presents a significant risk on organ survival.<sup>4</sup>
- Hepatic artery flows of less than 200 mL/min following orthotopic liver transplantation increase the risk of subsequent hepatic artery thrombosis six times.<sup>5</sup>

### Equipment Needed



HT364 Dual-channel Optima Flowmeter permits simultaneous measurements with two Flowprobes.



8-14 mm -AU Confidence Flowprobes® provide highly accurate measurements in vessels with fluctuating flows such as the portal vein. The Probes may be left in place for extended measurements and then easily removed via a ring attached to the pliable liner that cushions and protects the vessel.

### References

- <sup>1</sup>Henderson JM et al, "Hemodynamics During Liver Transplantation: The Interactions Between Cardiac Output and Portal Venous and Hepatic Arterial Flows," *Hepatology* 1992; 16(3): 715-718.
- <sup>2</sup>Henderson JM et al, "Volumetric and Functional Liver Blood Flow Are Both Increased in the Human Transplanted Liver," *J Hepatology* 1993; 17: 204-207.
- <sup>3</sup>Troisi R, de Hemptinne B, "Clinical Relevance of Adapting Portal Vein Flow in Living Donor Liver Transplantation in Adult Patients," *Liver Transplantation* 2004; 9(9): S36-S41.
- <sup>4</sup>Lin M et al, "Hepatic Artery Thrombosis and Intraoperative Hepatic Artery Flow Rates in Adult Orthotopic Liver Transplantation," *ANZ J Surg* 2002; 72: 798-800.
- <sup>5</sup>Pratschke S et al, "Arterial Blood Flow Predicts Graft Survival in Liver Transplant Patients," *Liver Transplantation* 2011; 17: 436-445.

Hashimoto K, Miller CM, Quintini C, Aucejo FN, Hirose K, Uso TD, Trenti L, Kelly DM, Winans CG, Vogt DP, Egtesad B, Fung JJ, "Is impaired hepatic arterial buffer response a risk factor for biliary anastomotic stricture in liver transplant recipients?" *Surgery* 2010; 148(3): 582-8. (9648AHM)

Kelly DM, Shiba H, Nakagawa S, Irefin S, Egtesad B, Quintini C, Aucejo F, Hashimoto K, Fung JJ, Miller C., "Hepatic blood flow plays an important role in ischemia-reperfusion injury." *Liver Transpl.* 2011 Dec;17(12):1448-56. (9647AHM)

### FLOWPROBE RECOMMENDATIONS

VESSEL	Probe Size (mm)	Probe Series
Hepatic artery	4 - 8	-FMV
Portal vein	8 - 14	-FMV, -AU
Common iliac a	8	-FMV, -FSB



4 mm and 6mm FMV Vascular Handle Flowprobes are recommended for hepatic arterial flow measurements.



8 to 14 mm FMV Vascular Handle Flowprobes are recommended for portal venous flow measurements.

# Adult Liver Donor Liver Transplantation

## Hepatic Artery & Portal Vein

### Protocol

### Living Donor

Measure right hepatic arterial and portal venous flow before hilar dissection.

Document measurements to serve as guide for expected flows in the recipient.

### Recipient

#### Recipient Hepatic Flow

Measure hepatic blood flow  
- following reperfusion  
- before biliary anastomosis  
- before wound closure

Compare with pre-transplant hepatic arterial flow

< 50 mL/min

Examine anastomosis for arterial thrombosis

Remeasure hepatic flow

Flow has increased

> 100 mL/min

Document flows and save waveforms for the operative record for post-op diagnostic consideration

#### Recipient Portal Flow

Measure portal blood flow  
- following reperfusion  
- after portal pressure measurement  
- before biliary anastomosis

Compare with pre-transplant portal venous flow

Flow increased up to 3 times pre-transplant portal flow

Flow increased > 3 times pre-transplant portal flow or >250 mL/min/110 gram graft weight

Reduced graft inflow by shunting portal flow away from liver<sup>1</sup>

Remeasure portal flow

Troisi R, de Hemptinne B, "Clinical Relevance of Adapting Portal Vein Flow in Living Donor Liver Transplantation in Adult Patients," Liver Transplantation 2004;9(9) Suppl 1 pp S36-S41. (6884AH)

## Hepatic/Portal References Cont.

Quintini C, Hirose K, Hashimoto K, Diago T, Aucejo F, Eghtesad B, Vogt D, Pierce G, Baker M, Kelly D, Miller CM, "Splenic artery steal syndrome" is a misnomer: the cause is portal hyperperfusion, not arterial siphon." *Liver Transpl.* 2008 Mar;14(3):374-9.(9649AHM)

Aucejo, FN, Hashimoto, K, Quintini, C, Kelly, D, Vogt, D, Winans, C, Eghtesad, B, Baker, M, Fung, J, Miller, C, "Triple-Phase Computed Tomography and Intraoperative Flow Measurements Improve the Management of Portosystemic Shunts during Liver Transplantation," *Liver Transplantation* 2008; 14: 96-99. (7606AH)

Wagener G, Gubitosa G, Renz J, Kinhabwala M, Brentjens T, Guarreram JV, Emond J, Lee HT, Landry D, "Vasopressin Decreases Portal Vein Pressure and Flow in the Native Liver during Liver Transplantation," *Liver Transplantation* 2008; 14: 1664-1670. (6673AH).

Aneman A., Eisenhofer G., Olbe L., Dalenback J., Nitescu P., Fandriks L., Friberg P., "Sympathetic Discharge to Mesenteric Organs and the Liver," *J Clin Invest* 1996; 97(5):1640-6. (1488AH)

Doi R, Inoue K., Kogire M., Sumi S., Takaori K., Suzuki T., Tobe T., "Simultaneous Measurement of Hepatic Arterial and Portal Venous Flows by Transit-time Ultrasonic Volume Flowmetry," *Surgery, Gynecology & Obstetrics* 1988;167(1):65-69. (26AH)

Doi, R., Inoue, K., Kogire, M., Sumi, S., Takaori, K., Suzuki, T., Tobe, T., "Study on Splanchnic Circulation: Measurement of the Liver Blood Flow," *Nippon Geka Gakkai Zasshi*, 1988;89(4):560-7. (109AH)

Figueras L., Llado L., Ramos E., Jaurrieta E., Rafecas A., Fabregat J., Torras J., Sabate A., Dalmau A, "Temporary portocaval shunt during liver transplantation with vena cava preservation. Results of a prospective randomized study," *Liver Transpl*, Vol. 7, No. 10, p. 904-11, 2001. (2201AH)

Panaro F, Bouyabrane H, Carabalona JP, Marchand JP, Jaber S, Navarro F, "Hepatic artery kinking during liver transplantation: survey and prospective intraoperative flow measurement," *J Gastrointest Surg.* 2012 Aug;16(8):1524-30. (9796AHE)

Jakab F, Rath, Z. Schmal F, Nagy P, Faller J, "Changes in Hepatic Hemodynamics Due to Primary Liver Tumours," *HPB Surgery* 1996; 9(4) 245-248. (854AH)

Jakab F, Rath Z, Schmal F, Nagy P, Faller J, "A New Method to Measure Portal Venous and Hepatic Arterial Blood Flow in Patients Intraoperatively" *HPB Surgery* 1996; 9(4) 238-243. (855AH)

Jakab F, Rath Z, Schmal F, Nagy P, Faller J, "The Afferent Circulation of the Liver in Patients with Primary Hepatocellular Carcinoma," *Hepatogastroenterology* 1995;42(4) 399-402. (704AH)

Jakab F, Rath Z, Schmal F, Nagy P, Faller J, "Blood Flow Measurement in Patients with Hepatocellular Carcinomas," *Acta Chir Hung* 1994;34(1-2): 87-94. (559AH)

Jakab F, Rath Z, Schmal F, Nagy P, Faller J, "The Interaction between Hepatic Arterial and Portal Venous Blood Flows; Simultaneous Measurement by Transit-Time Ultrasound Volume Flowmetry," *Hepatogastroenterology* 1995;42(1): 18-21. (560AH)

Jakab F, Rath Z, Schmal F, Nagy P, Faller J, "Intraoperative Estimation Bueno J, Escartin A, Baisells J, Margarit C, "Intraoperative Flow Measurement of Native Liver Allograft during Orthotopic Liver Transplantation in Children, *Transplant Proc.* 2007 39:7:2278-9. 7605AHM)

Rasmussen, A., Hjortrup, A., Kirkegaard, P., "Intraoperative Measurement of Graft Blood Flow - A Necessity in Liver Transplantation," *Transplant Int* 1997;10(1):774-77. (1761AHM)

Nanashima A, Pillay P, Crawford M, Nakasuji M, Verran DJ, Painter D, "Analysis of post-revascularization syndrome after orthotopic liver transplantation: the experience of an Australian liver transplantation center," *J Hepatobiliary Pancreat Surg*, Vol. 8, No. 6, p. 557-63, 2001. (2276AH)

Gontarczyk GW, Łągiewska B, Pacholczyk M, Trzebicki J, Jureczko L, Kołacz M, Kosieradzki M, Adamyński L, Wasiak D, Rowiński W, "Intraoperative blood flow measurements and liver allograft function: preliminary results," *Transplant Proc.* 2006 Jan-Feb;38(1):234-6.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

### AMERICAS

Transonic Systems Inc.  
34 Dutch Mill Rd  
Ithaca, NY 14850  
U.S.A.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

### EUROPE

Transonic Europe B.V.  
Business Park Stein 205  
6181 MB Elsloo  
The Netherlands  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

### ASIA/PACIFIC

Transonic Asia Inc.  
6F-3 No 5 Hangsiang Rd  
Dayuan, Taoyuan County  
33747 Taiwan, R.O.C.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

### JAPAN

Transonic Japan Inc.  
KS Bldg 201, 735-4 Kita-Akitsu  
Tokorozawa Saitama  
359-0038 Japan  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
info@transonic.jp