

# T400-Series Surgical Protocol

## Dog Hepatic Artery & Portal Vein: Chronic Blood Flow Measurement

### APPLICATION BASICS

Site: Hepatic Artery  
 Species: Dog  
 Weight: 25 kg  
 Duration: Chronic  
 Vessel Diameter: 2.5 mm

### PROBE

Size: 3 mm (side exit)  
 Reflector: L with sliding cover  
 Other: Silicone wrap  
 Catalog #: MC-3PSS-LS-WC100-CM4B-GC

**FLOWMETER** TS420 Perivascular Module

### APPLICATION BASICS

Site: Portal vein  
 Species: Dog  
 Weight: 25 kg  
 Duration: Chronic  
 Vessel Diameter: 12 mm

### PROBE

Size: 12 mm (side exit)  
 Reflector: U with Silicone Shield  
 Cable Length: 1 meter  
 Catalog #: MC-12PSS-USW-WC100-CRS10-GC

**FLOWMETER** TS420 Perivascular Module

## Application

The measurement of portal and hepatic blood flow has an important role in the study of hepatic physiology. One protocol was developed to study the hepatic extraction of metabolic substrates and also included the implantation of vascular occluders on the portal vein and the hepatic artery. Vascular access ports were also implanted in the hepatic vein, the hepatic artery and the portal vein. The concentration of any substrate may be determined from blood samples drawn from the vascular access ports. Since the total metabolic flux is the product of blood flow and the substrate concentration, total hepatic extraction may be determined. The relative contributions of the hepatic and portal vessels can be varied at will with the vascular occluders.

## Surgical Approach

Premedicate with 0.02 mg/kg atropine IM. Induce with 18 mg/kg thiamylal and maintain on halothane. With anesthetized dog in dorsal recumbency, make a midline skin incision from the xiphoid cartilage to the umbilicus. Continue the incision through the linea alba and the peritoneum to expose the lobes of the liver. Deflect the lobes of the liver cranially and identify the splanchnic vessels. Carefully dissect free a 2 cm segment of the portal vein and strip all fat from it for proper acoustical coupling. Slip the large U bracket around the vein, attach the body of the Probe and secure the screws. Rotation of the Probe around the vein may be necessary to align the screwdriver with each screw. Suture the cable to the perivascular connective tissue.

*(Continued on next side.)*

## Flow Ranges Observed



Fig. 1: Instantaneous flow ranged from 200 to 350 ml/min. This pulsatile flow example is from an acute experiment in an anesthetized dog. Note the periodic spikes from respiratory activity. Typical mean flow from chronic experiments ranged from 376 to 450 ml/min.

## Dog Hepatic Artery & Portal Vein: Chronic Blood Flow Measurement Cont.

### Surgical Approach cont.

Identify the gastroduodenal artery. It is the continuation of the hepatic artery after the last hepatic branch and runs adjacent to the bile duct. Ligate it to isolate the hepatic circulation from that of the stomach and pancreas.

Mobilize a 1 cm segment of the hepatic artery and strip away any fat. Pass the L bracket around the hepatic artery, close the slide and secure the screw. Wrap the silicone sheet around the Probe so that the parts of the bracket with the suture holes extend through the cutouts in the silicone sheet. Suture the ends of the silicone sheet as shown in Fig. 2. Place several 4-0 silk sutures between the each edge of the silicone wrap and perivascular connective tissue as shown in Fig. 3 Also place a single suture around the cable for strain relief.

Exit the Probe cables through a stab incision high on the abdominal wall. Make a skin incision between the shoulder blades and create a subcutaneous tunnel from the stab incision to the midscapular incision. Pull the cables through the subcutaneous tunnel. Close the body wall with 2-0 silk sutures in a simple interrupted pattern. Close the skin and the stab incisions with a subcuticular pattern.




Fig. 1: Flowprobes on hepatic artery & portal vein

Fig. 2: Silicone Wrap

Fig. 3: Suturing

### ACKNOWLEDGEMENT

Flow trace: Dr. Ronald Beck, University of South Carolina, School of Medicine, Dept. of Physiology, Columbia, SC 29208. ral Science, Uppsala, Sweden

## Dog Hepatic Artery & Portal Vein: Chronic Blood Flow Measurement Cont.

### REFERENCES

Eyer, S., Cerra, F., Konstantinides, F., Grim, E., "Splanchnic Hemodynamics and Extraction of Amino Acids Following a Meat Meal in the Awake Dog" FASEB Abstracts, 1990.

Beck R.R., Abel, F.L., "The Effects of ATP-Magnesium Chloride Infusion on Canine Endotoxin Shock", Circulatory Shock, 24:257, 1988.

Transonic Systems' Video "Measuring Hepatic Flow: Implanting the Transonic Flowprobe on the Hepatic Artery and the Portal Vein of the Dog" with Dr. Peter Alden.

Alden, P.B., Gonzalez, I., Morton, T., Chao, R.Y.N., Konstantinides, F.N., Cerra, B.F., "Hyperdynamic Sepsis Does Not Inhibit Hepatic Ketogenesis". Presented at the American Society of Critical Care Medicine annual meeting, Anaheim CA, 1987.

Kogire, M., Inoue, K., Doi, R., Sumi, S., Takaori, K., Suzuki, T., Tobe, T. "Effects of Intravenous Ethanol on Hepatic and Pancreatic Blood Flow in Dogs", Digest Dis Science, Vol. 33, No. 5, p. 592-97, 1988.

Nishiwaki, H., Ohira, M., Boku, T., Ishikawa, T., Nakagawa, H., Yata, K., Umeyama, K., "The Measurement of Hepatic Circulation before and after Orthotopic Liver Transplantation in Dogs — by Using Transit-Time Blood Flow Meter and H<sub>2</sub> Clearance Method", Nippon, Geka Gekka Zasshi, Vol. 90, No. 2, p. 243-9, 1989.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

#### AMERICAS

Transonic Systems Inc.  
34 Dutch Mill Rd  
Ithaca, NY 14850  
U.S.A.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

#### EUROPE

Transonic Europe B.V.  
Business Park Stein 205  
6181 MB Elsloo  
The Netherlands  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

#### ASIA/PACIFIC

Transonic Asia Inc.  
6F-3 No 5 Hangsiang Rd  
Dayuan, Taoyuan County  
33747 Taiwan, R.O.C.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

#### JAPAN

Transonic Japan Inc.  
KS Bldg 201, 735-4 Kita-Akitsu  
Tokorozawa Saitama  
359-0038 Japan  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
info@transonic.jp