

T400-Series Surgical Protocol

Fetal Sheep Umbilical Artery (Ventral Approach): Chronic Blood Flow Measurement

APPLICATION BASICS

Site:	Paired umbilical arteries
Species:	Fetal Sheep
State of Gestation:	109-120 days
Duration:	Chronic, 5 days
Vessel Diameter:	4 mm

PROBE


Size:	4 mm (side exit)
Reflector:	J
Connector:	10-pin
Cable Length:	60 cm
Catalog #:	MC-4PSS-JN-WC60-CRA10-GC

FLOWMETER

TS420 Perivascular Module

Flow Ranges Observed

Umbilical flow was 210 ml/min/kg fetal weight in two fetuses and 105 ml/min/kg in one fetus.



Application

Umbilical arterial blood flow is used extensively in pregnancy research. Some investigators use blood flow in combination with pressure to measure changes in vascular resistance induced by pharmaceutical agents. Others look for diurnal patterns associated with parturition.

Surgical Approach

Premedicate with 0.4 g glycopyrrolate IM, induce with 1 g ketamine IM and maintain anesthesia on 1.5% - 2% halothane.

Place anesthetized sheep in dorsal recumbency and make a ventral paramedian incision from the umbilicus to a point 2 cm cranial to the udder. The skin incision is made 1 cm off midline to avoid the median subcutaneous vein. Retract the skin and associated vascular structures and continue the incision through midline of the abdominal wall. Identify and exteriorize the uterine horn containing the fetus. Palpate the fetus to identify the orientation and make a transverse incision in the uterus to allow exteriorization of the fetus. Once the fetus is exposed make a midline skin incision just caudal to the umbilical cord to expose the paired arteries. Place the J bracket around the paired arteries and suture the Probe in place. Close the abdominal musculature and the fetal skin in separate layers with 2-0 simple continuous sutures. Secure the Probe cable to fetal skin with a 2-0 simple interrupted suture.

Close the uterus with a continuous Cushing pattern oversewn with a continuous Lembert, the Probe cable is exteriorized through this incision. Extend the Lembert pattern slightly to oversee the cable for 2 cm. Pass the Probe cable through the abdominal and close the body wall and skin routinely. Suture the cable to the skin of the ewe near the exit site.

Fetal Sheep Umbilical Artery (Ventral Approach): Chronic Blood Flow Measurement Cont.

ACKNOWLEDGEMENT

Kozo Akagi, Chikara Endo, Junya Saito, Masayuki Onodera, Masami Kozu, Shingo Tanigawara, Kunihiro Okamura, and Akira Yajima, Department of Obstetrics and Gynecology, Tohoku University School of Medicine, Sendai, 980, Japan.

Akira Sato, Department of Obstetrics and Gynecology, Fukushima Medical College, Fukushima, 960, Japan.

REFERENCES

Transonic Systems' 1989 FASEB Tutorial Video "Blood Volume Flow Measurement in the Fetus" with P Nathanielsz, J Dunn and R Wilkening

Akagi K, Endo C, Saito J, Onodera M, Tanigawara S, Okamura K, Yajima A, Sato A, "Ultrasonic Transit-Time Measurement of Blood Flow in the Animal Chronic Preparation Model," Jpn J Med Ultrasonics 1987; 14(2): 26-32.

Rudolph AM and Heymann MA 1980. Methods for studying the circulation of the fetus in utero. In: Monographs in Fetal Physiology: Animal Models in Fetal Medicine (I). Nathanielsz, P.W., editor. Pub. Perinatology Press, Ithaca, NY. pp. 1-58.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

AMERICAS

Transonic Systems Inc.
34 Dutch Mill Rd
Ithaca, NY 14850
U.S.A.
Tel: +1 607-257-5300
Fax: +1 607-257-7256
support@transonic.com

EUROPE

Transonic Europe B.V.
Business Park Stein 205
6181 MB Elsloo
The Netherlands
Tel: +31 43-407-7200
Fax: +31 43-407-7201
europe@transonic.com

ASIA/PACIFIC

Transonic Asia Inc.
6F-3 No 5 Hangsiang Rd
Dayuan, Taoyuan County
33747 Taiwan, R.O.C.
Tel: +886 3399-5806
Fax: +886 3399-5805
support@transonicasia.com

JAPAN

Transonic Japan Inc.
KS Bldg 201, 735-4 Kita-Akitsu
Tokorozawa Saitama
359-0038 Japan
Tel: +81 04-2946-8541
Fax: +81 04-2946-8542
info@transonic.jp