

Medical Note

Intraoperative Blood Flow Measurement during STA-M4/MCA Bypass Surgery for Moyamoya Revascularization

Introduction

One strategy to alleviate the symptoms of Moyamoya syndrome is the surgical creation of an arterial extracranial to intracranial (EC-IC) bypass from the superficial temporal artery (STA) to the M4 middle cerebral artery branch. The bypass is designed to augment flow in the intracranial territories (Fig. 1).

During surgery, the Charbel Micro-Flowprobe® is used to measure direct volume blood flow in the STA bypass and small target M4/MCA vessels. Intraoperative blood flow measurements confirm the quality of the anastomosis and assure that the target area is receiving sufficient blood from the bypass. Measurements also prompt revision if a technical error is suspected.

Flow Measurement Steps

Measure mean arterial pressure (MAP), end-tidal CO₂ and temperature. Record values on an Bypass Flow Record.

Pre-anastomosis: Intracranial Recipient Artery

1. Measure the diameter of the intracranial recipient artery (M4/MCA) and choose an appropriately sized Charbel Micro-Flowprobe® to measure recipient vessel flow.

PROBE SIZE	VESSEL RANGE, OUTER DIAMETER
1.5 mm	1.0 - 1.5 mm
2 mm	1.5 - 2.7 mm
3 mm	2.5 - 3.7 mm

2. Measure recipient vessel (M4/MCA) flow.
3. Record flow and flow direction on EC-IC Bypass Record.

Extracranial Donor Artery

4. Dissect the extracranial STA artery free, and skeletonize a segment for application of the Flowprobe.
5. Measure the diameter of the extracranial donor artery (STA) and choose the appropriately sized Flowprobe to measure STA baseline flow.

PROBE SIZE	VESSEL RANGE, OUTER DIAMETER
1.5 mm	1.0 - 1.5 mm
2 mm	1.5 - 2.7 mm
3 mm	2.5 - 3.7 mm
4 mm	3.3 - 4.7 mm
6 mm	4.4 - 6.6 mm




Fig. 1: Lateral view of an arterial flow augmentation STA-MCA EC-IC Bypass used to revascularize cranial territories for Moyamoya syndrome. The STA is exposed, cut and anastomosed to the MCA via an end-to-side anastomosis.

Post-anastomotic Flow Measurements

6. After construction of the STA-MCA bypass, measure post-anastomotic flows in the intracranial and extracranial arteries sequentially in the following order:
 - 1) distal M4/MCA (Fig. 3);
 - 2) proximal M4/MCA;
 - 3) distal STA;
 - 4) proximal STA.
7. If post-bypass flow in the recipient artery (sum of absolute values of distal and proximal M4/MCA flow) is not significantly above the pre-bypass flow, reexamine the anastomosis and the bypass for kinks or twists and redo, if necessary. Apply a vasodilator (papaverine) when there has been some vasospasm due to manipulation of the vessel and/or flow measurements seem to be low or absent.
8. Record flow rates and flow directions, MAP, end-tidal CO₂, and occlusion time on the EC-IC Bypass Record.

Intraoperative Blood Flow Measurement during STA-MCA Bypass Surgery for Moyamoya Revascularization Cont.




Fig. 2: Photo shows the M4/MCA site just before the Flowprobe is slipped around the vessel to measure baseline M4 flow before anastomosing the bypass to the vessel. The blue background is placed to help visibility during sewing the anastomosis and as the Flowprobe is applied to the vessel.





Fig. 3: Measuring blood flow in recipient M4/MCA artery after anastomosis to STA bypass.

Protocol: Flow Measurement during EC-IC Bypass Revascularization for Moyamoya Syndrome



References

Lee M, et al, "Intraoperative blood flow analysis of direct revascularization in patients with moyamoya disease," J Cereb Blood Flow & Metab 2011;31(1):262-74. (7969AHAH)

Lee M, et al, "Quantitative hemodynamic studies in moyamoya disease: a review," Neurosurg Focus 2009; 26(4):E5. (6762AH)

Guzman R, et al, "Clinical Outcome after 450 Revascularization Procedures for Moyamoya Disease," J Neurosurg 2009; May 22.

Veeravagu A, Guzman R, Patil CG, Hou LC, Lee M, Steinberg GK "Moyamoya disease in pediatric patients: outcomes of neurosurgical interventions," Neurosurg Focus 2008;24(2):E16.

Chang, SD, Steinberg, GK, "Surgical Management of Moyamoya Disease," <http://www.moyamoya.com/journals/moyamoya.html>

Golby, AJ et al, "Direct and Combined Revascularization in Pediatric Moyamoya Disease," 1999;45:50-60.

Kawaguchi, S et al, "Effect of Direct Arterial Bypass on the Prevention of Future Stroke in Patients with the Hemorrhagic Variety of Moyamoya Disease." J Neurosurg, 2000;93:397-401. (7248AH)

Zipfel, GJ, et al, "Moyamoya Disease in Adults: The Role of Cerebral Revascularization." Skull Base, 2005; 1:27-41.