

# T400-Series Surgical Protocol

## Rat Mesenteric Branches: Acute Blood Flow Measurement

### APPLICATION BASICS

Site:	Mesenteric branches: jejunal & ileal arcades
Species:	Rat
Body Weight:	400 grams
Duration:	Acute
Vessel Diameter:	270 $\mu$ m
<b>PROBE</b>	
Size:	0.5 mm
Reflector:	J
Connector:	CRA10: 10-pin
Cable Length:	60 cm
Catalog #:	MA-0.5PSB
<b>FLOWMETER</b>	
	TS420 Perivascular Module

### Application

The combination of this model and the new microcirculation transit-time Probe, allowed us to make the first measurements of true flows in the microcirculation. These were until now estimated from diameters and velocity measurements; this method lacked precision and prevented any appreciation of the instantaneous dynamic of flow. This allowed us to examine the instantaneous dynamic and to observe the cyclic variation of flow caused by the periodic contractions of the bowel.

### Surgical Approach

Note: This protocol is for mesenteric branches, jejunal and ileal arcades. The superior mesenteric artery (0.7 - 0.9 mm diameter) is larger than the branches and requires 1 mm Probes for full ultrasonic illumination.

Male Sprague-Dawley rats weighing approximately 400 gm were anesthetized with 0.1 ml / 100 gm pentobarital (6%) i.p. Anesthesia was maintained by additional i.p. injections of diluted pentobarital. The right jugular vein was cannulated with polyethylene tubing for the administration of drugs.

After a small abdominal incision, a section of the ileum was pulled out and spread over the transparent stage. Fat and connective tissues surrounding the mesenteric arteries (internal artery diameters were approximately 270  $\mu$ m) were carefully removed under a dissecting microscope.

### Flow Ranges Observed




Fig. 1: Mean flow over 90 seconds was  $0.494 \pm 0.276$  and  $0.527 \pm 0.233$  ml/min, respectively. After treatment with isoprenaline and superfusion of the preparations with papaverine, mean flow increased to  $0.666 \pm 0.379$  ml/min.




Fig. 2: Transit-time ultrasound Flowprobe on mesenteric (jejunal) branch (center); Laser Doppler Probe on mucosa on the right.

## Rat Mesenteric Branches: Acute Blood Flow Measurement Cont.

### Surgical Approach cont.

The preparation was mounted under a biological binocular microscope (Leitz) connected to a color video recording system. The surface of the mesenteric artery was covered with a Saran film. Warm Krebs' solution was superfused on the artery at a rate of 2 ml/min. The microscope magnification was 145 fold. The stainless steel handle of the Probe was connected to a micromanipulator in order to allow positioning of the Probe around the vessel near the objective. The vessel was positioned in the lumen of the Probe, making sure the reflector bracket was not tugging on the vessel wall to reduce the flow. Zero flow reading was obtained by transient clamping of the arteriole under investigation.

Basal flow measurements were done at a few minute intervals, the mean flows over 90 seconds were  $0.494 \pm 0.276$  ml/min and  $0.527 \pm 0.233$  ml/min respectively. The reproducibility seems excellent despite the elevated standard deviation due to the physiological periodic reduction of flow provoked by intestinal peristalsis.

After treatment of the animals with isoprenaline and superfusion of the preparations with papaverine, the mean measured flow increased to  $0.666 \pm 0.379$  ml/min.

### REFERENCES

Marque S et al, "APC Inhibits FMLP Induced Leukocytes Sticking in the Micro-circulation," *FASEB J* 1993; 7(3.1): A344.

#### Superior Mesenteric Blood Flow

Myers SI et al, "Oxygen Free Radicals Regulate Splanchnic Nitric Oxide Synthesis and Blood Flow," *Cardiovascular Surgery* 1995; 3(2): 207-10.

Wang J F et al, "The Roles of Nitric Oxide and Hydrogen Peroxide Production in Lipopolysaccharide-Induced Intestinal Damage," *Shock* 1995; 2(3): 185-191.

Turnage R H et al, "Neutrophil Regulation of Splanchnic Blood Flow after Hemorrhagic Shock," *Ann Surg* 1995; 222(1) 66-72. 575A)

Pofahl, W.F. et al, "Small Intestinal Microcirculatory Effects of Octreotide," *J Surg Res* 1994; 56:345-350. (397A)

Kost C K et al, "Vascular Reactivity to Angiotensin II is Selectively Enhanced in the Kidneys of Spontaneously Hypertensive Rats," *J Pharmacol & Expe Therap* 1994; 269(1) 82-8. (426A)

Jackson EK, Herzer WA, "Angiotensin II/Prostaglandin I2 Interactions in Spontaneously Hypertensive Rats," *Hypertension* 1993, 22: 688-98. (344A)

Turnage RH et al, "Splanchnic PG12 Release and "No Reflow" Following Intestinal Perfusion," *J Surg Res* 1995; 222(1): 66-72.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

#### AMERICAS

Transonic Systems Inc.  
34 Dutch Mill Rd  
Ithaca, NY 14850  
U.S.A.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

#### EUROPE

Transonic Europe B.V.  
Business Park Stein 205  
6181 MB Elsloo  
The Netherlands  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

#### ASIA/PACIFIC

Transonic Asia Inc.  
6F-3 No 5 Hangsiang Rd  
Dayuan, Taoyuan County  
33747 Taiwan, R.O.C.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

#### JAPAN

Transonic Japan Inc.  
KS Bldg 201, 735-4 Kita-Akitsu  
Tokorozawa Saitama  
359-0038 Japan  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
info@transonic.jp