

T400-Series Surgical Protocol

Rat Renal Artery: Acute Blood Flow Measurement

APPLICATION BASICS

Site:	Renal artery
Species:	Rat
Body Weight:	< 300 grams
Duration:	Acute
Vessel Diameter:	0.7 - 0.8 mm

PROBE

Size:	1 mm
Reflector:	V / JS
Connector:	CRA10: 10-pin
Cable Length:	60 cm
Catalog #:	MA-0.7V MA-1PRB

FLOWMETER	TS420 Perivascular Module
------------------	---------------------------

Application

The measurement of renal blood flow has an important role in research on hemodynamics, electrolyte regulation and pregnancy induced hypertension. Flow-pressure relationships are essential in defining renal autoregulation. Other studies have focused on diuretics, cardiovascular drugs, and nephrotoxic agents. While average renal flow may also be obtained from the renal vein, the pulsatile waveform of the renal artery provides additional information and visual confirmation of a measurement with a properly functioning Flowprobe.

Surgical Approach

Anesthetize the rat with Inactin anesthesia (100 mg/kg IP). If Inactin is not available, sodium pentobarbital (60 mg/kg IP) may be used instead. Note that pentobarbital anesthesia is less stable than Inactin and may significantly depress flow when compared to that of the conscious animal. The use of a heating pad or hot water bottle is recommended as hypothermia also reduces flow. In long procedures, fluid infusion (0.9% NaCl @ 1 ml/hr) through a femoral catheter is also recommended.

Place rat in dorsal recumbency and make a ventral midline abdominal skin incision. Extend the abdominal incision through the linea alba into the abdominal cavity. Deflect the intestines to the rat's right to expose the left kidney. To expose the right kidney, deflect the intestines to left. Identify the large renal

(Continued on next side.)

Flow Ranges Observed




Fig. 1: Renal flow in the rat varies greatly with the anesthetic used and the plane of anesthesia. In general, protocols with ketamine will show higher flows than those with pentobarbital. Hypothermia is also a common cause of lower than expected flow measurements.




Fig. 2: Flowprobe on renal artery.

Rat Renal Artery: Acute Blood Flow Measurement Cont.

Surgical Approach cont.

vein; the renal artery is much smaller, cranial and deep to the vein. Carefully locate and dissect out the renal artery. Remove adjacent fat for proper acoustical coupling. Place the Probe around the artery (Fig. 2). Manually position the artery so that the vessel lies within the lumen of the Probe. A micromanipulator may be used to stabilize the Flowprobe position on the artery.

Remove the plunger of a 30 cc syringe and load the syringe with Surgilube acoustic gel, taking care to prevent the formation of air bubbles. Place a flexible catheter on the tip of the syringe. Insert the flexible catheter through the Probe's acoustic window adjacent to the artery and deposit the gel while withdrawing the syringe. The lubricating gel acts as an acoustical couplant and must replace all air space. Check the signal bar indicator on the Meter to verify that signal amplitude is about 1 Volt. A low signal or an acoustic error can usually be traced to an insufficient amount of lubricating gel or an air bubble.

ACKNOWLEDGEMENTS

Dr. Brian Murray, Department of Medicine, Nephrology Division, State University of New York at Buffalo, Buffalo, NY

Dr. Thomas L. Smith, Department of Orthopedic Surgery, Wake Forest University School of Medicine, Winston-Salem, NC.

REFERENCES

Welch WJ et al, "Validation of Miniature Ultrasonic Transit Time Flow Probes for Measurement of Renal Blood Flow in Rats," AJP 1995; 268,(1)Pt 2, F175-178.

Welch WJ et al, "The Ultrasonic Transit-time Blood Flow Meter: An Accurate Method for Measuring Renal Blood Flow in the Rat," JASN 1993; 3: 528.

Griffin KA et al, "Role of Endothelium-Derived Nitric Oxide in Hemodynamic Adaptations after Graded Renal Mass Reduciton," AJP 1993; 264: R1254-R1259.

Hatton DC et al, "Impact of Stress-Induced Hypertension on Renal Blood Flow and Vascular Resistance," JASN 1993; 4(3): 512.(336A)

Holycross BJ et al, "Adenosine Receptor Blockade Attenuates Angiotensin II Induced Decreases in Renal Blood Flow in Vivo," FASEB J 3(3): Renal Hemodynamics A109.

Hoffman A et al, "Endothelin Induces An Initial Increase in Cardiac Output Associated with Selective Vasodilation in Rats," Life Sciences 1989; 45(3): 249-255.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

AMERICAS

Transonic Systems Inc.
34 Dutch Mill Rd
Ithaca, NY 14850
U.S.A.
Tel: +1 607-257-5300
Fax: +1 607-257-7256
support@transonic.com

EUROPE

Transonic Europe B.V.
Business Park Stein 205
6181 MB Elsloo
The Netherlands
Tel: +31 43-407-7200
Fax: +31 43-407-7201
europe@transonic.com

ASIA/PACIFIC

Transonic Asia Inc.
6F-3 No 5 Hangsiang Rd
Dayuan, Taoyuan County
33747 Taiwan, R.O.C.
Tel: +886 3399-5806
Fax: +886 3399-5805
support@transonicasia.com

JAPAN

Transonic Japan Inc.
KS Bldg 201, 735-4 Kita-Akitsu
Tokorozawa Saitama
359-0038 Japan
Tel: +81 04-2946-8541
Fax: +81 04-2946-8542
info@transonic.jp