

# Scisense PV Technical Note

## Understanding Preload

Preload is known as the load imposed on the ventricle at the end of diastole. At a cellular level, preload is defined as the maximum degree of myocardial fiber stretch or tension before ventricular contraction, determined by the mean sarcomere length at the end of diastole. Since sarcomere length cannot be determined in the intact heart, other indices of preload such as ventricular end diastolic volume (EDV) or pressure (EDP) are used. In general, EDV offers a better estimation of preload than EDP.


The relationship between the changes in preload and stroke volume depends on the morphology and Frank-Starling curve, which are determined by the contractile capacity of the heart and the ventricular afterload. Increasing preload increases stroke volume by a non-linear relationship. Cardiac preload can also be described as the passive filling properties of ventricles.

### FACTORS CAUSING PRELOAD (EDV) INCREASE

- Increased ventricular compliance
- Venoconstriction
- Increased venous return (skeletal muscle activity and respiratory activity during physical activity or position and gravity)
- Decreased heart rate (increased filling time)
- Neuro-endocrine stimulation of venous tone increases (excitement)
- Increased blood volume (e.g. post-transfusion)


### FACTORS CAUSING PRELOAD (EDV) DECREASE

- Compliance of ventricle decreases or stiffness increases (multiple muscle diseases leading to hypertrophy, post-tissue graft implantation)
- Increased heart rate (reduced filling time)
- Venodilatation (peripheral or central) causing blood to pool in legs, abdomen, liver etc.
- Atrial arrhythmias (impaired contraction)



Different PV loops are obtained with different preloads, modeled by constant contractility (ESPVR and EDPVR boundaries) and afterload.

### Dilated



During Dilated Cardiomyopathy Frank-Starling curves shifts down and to the right due to a decrease in contractility. During this dysfunction, stroke volume is reduced and preload increases as a consequence. The final result is a more compliant LV chamber, as seen by the increase in ESV and decreased slope/rightward shift of ESPVR.

## Understanding Preload Cont.

Preload is described by end diastolic volume and can be measured by a variety of methods.

### THERMODILUTION

Thermodilution can be used to obtain a global EDV index to evaluate biventricular preload. Alternatively, right ventricle EDV can be obtained using a pulmonary artery catheter with a rapid response thermistor in the right ventricle.


### ECHOCARDIOGRAPHY

Echocardiography (both transthoracic and transesophageal) gives reliable measures of end diastolic surface areas. End-diastolic volume can then be calculated using Simpson's method or a similar approach. Long axis Doppler echocardiography can measure increases in stroke volume (SV) during the period of increased preload.


### PRESSURE-VOLUME LOOPS

Pressure-volume loops give a direct measurement of end diastolic volume as derived by admittance. PV loops have the advantage of also providing information on heart contractility based on the end systolic and end diastolic pressure volume relationships (ESPVR and EDPVR). Additionally, PV loops provide stroke volume and end diastolic pressure measurements for determining the Frank-Starling curve.

### Hypertrophic



### Restrictive



**Both left ventricular hypertrophy and restrictive cardiomyopathy exhibit a reduced stroke volume and thus preload. The LV becomes less compliant as seen by the leftward shift in EDPVR.**



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

#### AMERICAS

Transonic Systems Inc.  
34 Dutch Mill Rd  
Ithaca, NY 14850  
U.S.A.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

#### EUROPE

Transonic Europe B.V.  
Business Park Stein 205  
6181 MB Elsloo  
The Netherlands  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

#### ASIA/PACIFIC

Transonic Asia Inc.  
6F-3 No 5 Hangsiang Rd  
Dayuan, Taoyuan County  
33747 Taiwan, R.O.C.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transonicasia.com

#### JAPAN

Transonic Japan Inc.  
KS Bldg 201, 735-4 Kita-Akitsu  
Tokorozawa Saitama  
359-0038 Japan  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
info@transonic.jp