

Comprehensive genomic profiling without compromises

The enhanced Ion Torrent™ Oncomine™ Comprehensive Assay Plus

- **Single-gene biomarkers**—detect all types of single-gene variants, such as single-nucleotide variants (SNVs), insertions and deletions (indels), novel and known fusions, splice variants, and copy number variants (CNVs), including both copy number gains and losses
- **Multiple-gene biomarkers**—study potential responses to immunotherapies with measurement of tumor mutational burden (TMB) and predisposition to genetic hypermutability by comparing microsatellite instability (MSI) regions, and analyze mutational signatures for insights into etiological factors in tumorigenesis
- **Homologous recombination repair deficiency (HRD) research**—detect both gene-level and sample-level loss of heterozygosity (LOH) to assess genomic instability and mutations in 46 key genes in the homologous recombination repair (HRR) pathway
- **Low input requirements**—formalin-fixed, paraffin-embedded (FFPE) sample inputs of 20 ng DNA or RNA are sufficient to profile more than 500 genes, helping ensure more samples can be analyzed
- **High testing success**—high sequencing success rates (up to 95%), combined with low QNS (quantity not sufficient) results, help ensure more samples are successfully tested
- **Bioinformatics solution**—streamlined bioinformatics analysis pipeline is optimized for this assay and packaged in a user-friendly experience with the Ion Torrent™ Oncomine™ Reporter, which gives fully annotated results
- **Highly automated workflow**—hands-on time of ~1 hour supports lab efficiency and helps to reduce possible errors due to handling

End-to-end solution

Automated library and template preparation on the Ion Chef™ System

Sequencing on the Ion GeneStudio™ S5 System

End-to-end bioinformatics solution

Oncomine Reporter generates fully annotated results



40 min of hands-on time



20 min of hands-on time



Automated bioinformatics pipeline (7 hr of data processing)



<5 min for report generation

ion torrent

FusionSync technology

Key considerations for optimal fusion detection

- Fusion detection from low-input samples
- Detection of low level of fusion transcripts
- Ability to detect novel fusions for driver genes

With Ion Torrent™ FusionSync™ technology, the Oncomine Comprehensive Assay Plus now covers >1,300 isoforms across 49 fusion drivers, enabling highly sensitive and robust detection of known fusions and novel combinations of known fusion partners at low levels of a tumor-specific fusion transcript in a background of normal RNA and with minimal

input down to 20 ng. Simultaneously, the exon-tiling imbalance approach is now enabling this assay to detect novel fusions with key fusion driver genes, such as *ALK*, *FGFR2*, *NTRK1*, *NTRK2*, *NTRK3*, and *RET* (Figure 1). For each driver gene in which a fusion is detected, the software also predicts with high confidence the position of the fusion breakpoint relative to the kinase domain. This is critical, as an intact kinase domain is essential for the pathogenicity of a fusion event.




Figure 1. Detecting a novel fusion via an exon-tiling imbalance approach. Source: Internal R&D data.

- Software measures the intragenic 3' to 5' expression ratio for each gene and compares the ratio to the baseline (normal sample)
- Genes that do not undergo a fusion event are expected to have a 3' to 5' expression ratio similar to the baseline, i.e., an imbalance score close to zero
- Genes that undergo a fusion event typically have a 3' to 5' expression ratio greater than the baseline, i.e., an imbalance score much greater than zero

HRD research

HRD is becoming an important biomarker in precision oncology clinical research. Under normal conditions, errors during homologous recombination are repaired in the HRR pathway (Figure 2). Errors in the HRR pathway, such as loss-of-function or deleterious mutations in the associated genes, lead to higher levels of genomic instability. The Oncomine Comprehensive Assay Plus is a single assay for comprehensive detection of HRR gene mutations, as well as detection of genomic instability by assessing LOH at both the gene level and sample level, and identifying mutational signatures.




Figure 2. HRR pathway. Non-grey genes are covered in the Oncomine Comprehensive Assay Plus. Purple genes were included in clinical trials with prostate cancer clinical research samples.

Mutation detection of HRR pathway genes

The significant role of HRR genes in maintaining genome stability and tumor suppression has been studied extensively, especially in the *BRCA1* and *BRCA2* genes. In recent years, it has been demonstrated that alterations in a *BRCA*ness pathway, including HRR genes, may increase the risk of developing tumors. The status of HRR genes is now considered a new potential biomarker for precision oncology. Figure 3A shows the HRR genes covered in this assay.

Genomic instability measurement

The Oncomine Comprehensive Assay Plus measures genomic instability with both gene-level and sample-level LOH with high accuracy. Figure 3 demonstrates the LOH assessment at both sample level and gene level with the Oncomine assay compared with the Applied Biosystems™ OncoScan™ CNV Assay as an orthogonal test using the same FFPE samples.

A
46 HRR Pathway genes are covered by the Oncomine Comprehensive Assay Plus

<i>ABRAXAS1</i>	<i>ATM</i>	<i>ATR</i>	<i>BAP1</i>	<i>BARD1</i>	<i>BLM</i>	<i>BRCA1</i>	<i>BRCA2</i>	<i>BRIP1</i>	<i>CDK12</i>	<i>CHEK1</i>	<i>CHEK2</i>
<i>FANCA</i>	<i>FANCC</i>	<i>FANCD2</i>	<i>FANCE</i>	<i>FANCF</i>	<i>FANCG</i>	<i>FANCI</i>	<i>FANCL</i>	<i>FANCM</i>	<i>MRE11</i>	<i>NBN</i>	<i>PALB2</i>
<i>PARP1</i>	<i>PARP2</i>	<i>PARP3</i>	<i>POLD1</i>	<i>POLE</i>	<i>PPP2R2A</i>	<i>PTEN</i>	<i>RAD50</i>	<i>RAD51</i>	<i>RAD51B</i>	<i>RAD51C</i>	<i>RAD51D</i>
<i>RAD52</i>	<i>RAD54L</i>	<i>RNASEH2A</i>	<i>RNASEH2B</i>	<i>RNASEH2C</i>	<i>RPA1</i>	<i>SLX4</i>	<i>TP53</i>	<i>XRCC2</i>	<i>XRCC3</i>		





Figure 3. HRR pathway research with Oncomine Comprehensive Assay Plus. (A) The 46 HRR pathway genes covered in the Oncomine Comprehensive Assay Plus. (B) Oncomine Comprehensive Assay Plus sample-level percent LOH estimates (y-axis) correlate favorably with the orthogonal test (x-axis), the OncoScan assay, on the same FFPE samples. Test sample set consists of FFPE samples from various solid-tumor tissue types. Pearson correlation (R) is shown as the measure of association. (C) Gene-level LOH accuracy comparing 20 genes in the HRR pathway (concordance defined as proportion of these 20 genes that have LOH in both OncoScan assay and Oncomine Comprehensive Assay Plus, with 90% mean accuracy across clonal samples). Source: Internal R&D data.

New: Bioinformatics solution enabling visualization of mutational signatures

Mutational signatures are an important tool for precision oncology research, providing insights into the biological mechanisms involved in carcinogenesis (e.g., UV damage, deficiency in DNA repair).


The Oncomine Comprehensive Assay Plus provides you with a comprehensive genomic profile, and as part of the streamlined bioinformatics workflow, the mutational signature plot is automatically generated and does not require additional analysis with third-party software.

Oncomine Comprehensive Assay Plus performance across single-biomarker variant types

Using commercially available reference controls and clinical research FFPE samples, assay sensitivity and specificity ranged from 93% to 100%, with averages of 97.0% sensitivity and 98.3% specificity across all variants, with CNV gain and CNV loss demonstrating exceptional 100% specificity.

Ordering information

Description	Quantity	Cat. No.
Oncomine Comprehensive Assay Plus, automated library preparation	32 rxn	A49667
Oncomine Comprehensive Assay Plus RNA, automated library preparation	32 rxn	A49671
Oncomine Comprehensive Assay Plus, manual library preparation	24 rxn	A48577
Oncomine Comprehensive Assay Plus RNA, manual library preparation	24 rxn	A48578

 Find more information about Ion Torrent™ Oncomine™ NGS solutions at oncomine.com/

ion torrent