


Your Guide to Life365's Online Portal



CARE COORDINATOR: NOTIFICATIONS WORKFLOW
FOR BREACHED THRESHOLD PARAMETERS

Coordinator Portal Core Functions

- 1 Review Member Populations Health Data
- 2 Establish Vital Thresholds for Members
- 3 Review Notifications for Breached Threshold Parameters
- 4 Assign, Remove, & Reassign Smart Devices
- 5 Graphing & Trending Functions
- 6 Create Reminders for Members
- 7 Review Survey Responses

Contents

Member List Dashboard	1
Notifications Tab	2
Notifications Workflow	3
Proceed Prompt	3
Member Vitals Page	4
Review Breached Threshold Parameter	5
Record Comments	6
Quick-view Recent Comments	7
Flag Breached Measurement	8

Member List Dashboard

The screenshot shows the Life365 Member List Dashboard. At the top left is the Life365 logo. To its right are navigation links for 'Support', 'English', and 'Life365 Coordinator'. Below the logo are two tabs: 'Vitals' and 'Notifications 12'. Callout 'B' points to the 'Notifications 12' tab. Below the tabs is a search bar with 'Search Users', 'Go', and 'Clear' buttons. Callout 'A' points to the 'Members List' section header. Below this is a table with columns 'Member' and 'Last Measurement'. The first row shows 'Jen Smith' with various vital signs: BP (112/75), ECG (N/A), Glc (106 mg/dL), HR (N/A), SpO2 (97%), Respiration Rate (N/A), SP (N/A), Temp (97.7 °F), and Wgt (145.3 lbs).

Member	Last Measurement ▼					
Jen Smith Coordinator: Emily Fellows Group: Life365, Inc. Site: Tempe IMEI: 123456789012345	Blood Pressure, 02/05/2020 09:42					
	BP 112/75 (85) 05 Feb 09:42	ECG N/A	Glc 106 mg/dL 15 Mar 10:09	HR N/A	SpO2 97% (70) 16 Jan 09:04	Respiration Rate N/A
	SP N/A	Temp 97.7 °F 06 Dec 08:50	Wgt 145.3 lbs 05 Feb 09:41			

Upon Sign-in, you will be directed to your Members List Dashboard (A).

Review the number of Notifications for breached threshold parameters within the Tab Notifications (B). To proceed with the Notifications workflow, Select the Notifications Tab (B).

Notifications Tab

Select Member Name to proceed with workflow.

Mouse-over the Member to reveal Member, measurement type breach, and date & time received.


Member	Received at	Details
Jen Smith Coordinator: Life365 Coordinator Group: Life365.Health Site: Scottsdale, AZ	Today 14:53	Oxygen Saturation Breach: 98 % (96 bpm) Oxygen Saturation Thresholds Oxygen Saturation 90 – 91 % Pulse 60 – 90 bpm
Michaela James	February 07, 2020 15:15	Oxygen Saturation Breach: 95 % (88 bpm)
Roger McDonald	February 07, 2020 15:02	Oxygen Saturation Breach: 97 % (82 bpm)

Displays Date & Time notification was received.

Display breached threshold details including:

- Measurement Type & Breach
- Threshold Parameter for Measurement Type

Notifications Workflow: Proceed Prompt



Selecting Member name will display prompt.

Select **Continue** to continue to Member vitals.

Check **Do not show again** to remove future prompts.

Notifications		
Member	Received at	Details
Jen Smith Coordinator: Life365 Coordinator Group: Life365.Health Site: Scottsdale, AZ	Today 14:53	Oxygen Saturation Breach: 98 % (96 bpm) Oxygen Saturation Thresholds Oxygen Saturation 90 – 91 % Pulse 60 – 90 bpm
Michaela James	February 07, 2020 15:15	Oxygen Saturation Breach: 95 % (88 bpm)
Roger McDonald	February 07, 2020 15:02	Oxygen Saturation Breach: 97 % (82 bpm)

Notifications Workflow: Member Vitals Page

Member Detail

Jen Smith, 22 yrs

- Profile
- Management Details
- Notes
- Vitals** (A)
- History
- Export
- Notifications
- Display Options

Blood Pressure
112/72 mmHg
78 bpm
11/19/2019 09:05

Oxygen Saturation (B)
98%
85 bpm
11/26/2019 10:16

ECG
No measurements at the moment.

(C)

You will be directed to the Members Vitals page (A).

Select the Measurement Type with the breached threshold parameter (B).

Note that breached measurement types will display a colored circle icon reflecting the Moderate or High breach (C).

Notifications Workflow: Review Breached Threshold Parameter

Date	Systolic, mmHg	Diastolic, mmHg	Pulse, bpm	Status
02/18/2020 10:04	96	58	76	●
02/18/2020 10:02	109	73	70	●
02/18/2020 09:58	109	73	71	N/A
02/18/2020 09:41	123	78	80	N/A
02/17/2020 14:11	122	70	72	N/A
02/17/2020 14:10	126	71	66	N/A

Select the vitals reading to review:

- Member Information
- Bluetooth Device Information
- Historical Comments
- Create Comments
- Flag Members Vitals Reading

Quick-view the Status of vital readings within the Status column.

- Red reflects High Breach
- Orange reflects Moderate Breach

Notifications Workflow: Record Comments

The screenshot displays a patient record for Jen Smith, 22 years old. The interface includes a sidebar with navigation options: History, Timeline, Blood Pressure, ECG, Blood Glucose, Heart Rate, Oxygen Saturation, Spirometry, Temperature, Weight, and Respiration Rate. The main content area is titled "Oxygen Measurement Details" and shows a reading of "97% 84 bpm". Below this is a section for "Member Details" and "Device Details". The Member Details section includes fields for Phone, Email Address, Preferred Language, and Site Location. The Device Details section includes fields for Name, Serial, Address, Manufacturer, Model, Firmware, and Hardware. At the bottom of the main content area, there is a "Add new comment" section with a text input field and a "Save" button. Below the comment input field, there is a "Yes" button and a "No" button. At the very bottom, there is a list of historical comments, with the first comment by Emily Fellows, updated 21 days ago, stating "Called Jen and received no response. Left voicemail for Jen to reach out when possible. Vital reading is Flagged for follow-up."

View vital reading.

Quick-view Member Details and Device Details.

Add a new comment in the provided space.

Select Save only once comment has been included.

Review historical comments from Coordinators below.

Notifications Workflow: Quick-View Recent Comments

The screenshot displays a 'Blood Pressure' history table with columns for Date, Systolic, mmHg, Diastolic, mmHg, Pulse, bpm, and Status. A 'Details' button is highlighted in an orange box. A callout window shows a detailed view of a comment, including the coordinator's name 'Emily Fellows', the date '01/28/2020 14:06', and the comment text: 'Called Jen and received no response. Left voicemail for Jen to reach out when possible. Vital reading is Flagged for follow-up.'

Date	Systolic, mmHg	Diastolic, mmHg	Pulse, bpm	Status	
02/18/2020 10:04	96	58	76	●	Details
02/18/2020 10:02	109	73	70	●	Details
02/18/2020 09:58					
02/18/2020 09:41					
02/17/2020 14:11					

Date	Systolic, mmHg	Diastolic, mmHg	Pulse, bpm	Status	
02/18/2020 10:04	96	58	76	●	Details
Note Created By: Emily Fellows Called Jen and received no response. Left voicemail for Jen to reach out when possible. Vital reading is Flagged for follow-up.					
01/28/2020 14:06					
02/18/2020 10:02	109	73	70	●	Details
02/18/2020 09:58	109	73	71	N/A	Details
02/18/2020 09:41	123	78	80	N/A	Details
02/17/2020 14:11	122	70	72	N/A	Details
02/17/2020 14:10	126	71	66	N/A	Details

Quick-view most recent comment for vital reading by Selecting the Details button.

Recent comment will display along with:

Name of Coordinator who created comment
Date & Time comment was created

Notifications Workflow: Flag Breached Measurements

Jen Smith, 22 yrs

History

- Timeline
- Blood Pressure
- ECG
- Blood Glucose
- Heart Rate
- Oxygen Saturation
- Spirometry
- Temperature
- Weight
- Respiration Rate

Oxygen Measurement Details

97% 84 bpm

Member Details	Device Details
Phone: (555) 123-4567	Name:
Email Address: jsmith@life365.health	Serial:
Preferred Language: Spanish	Address:
Site Location: Arizona	Manufacturer:
	Model:
	Firmware:
	Hardware:

Add new comment

? Yes No

Emily Fellows
Updated at 21 days ago
Called Jen and received no response. Left voicemail for Jen to reach out when possible. Vital reading is Flagged for follow-up.

Flag a Members breached vital reading by enabling Yes for flagging.

Selecting Yes will return the breached vital reading to the Notifications Tab.


NOTE: A comment must be created before saving a Flag change.

Enable **No** to save comment only.

Select **Save** when updated.

Notifications Workflow: Flag Breached Measurements

Notifications		
Member	Received at	Details
Jen Smith Coordinator: Life365 Coordinator Group: Life365.Health Site: Scottsdale, AZ	Today 14:53	Oxygen Saturation Breach: 98 % (96 bpm) Oxygen Saturation Thresholds Oxygen Saturation 90 – 91 % Pulse 60 – 90 bpm
Michaela James	February 07, 2020 15:15	Oxygen Saturation Breach: 95 % (88 bpm)
Roger McDonald	February 07, 2020 15:02	Oxygen Saturation Breach: 97 % (82 bpm)



In the Notifications Tab, mouse-over the Flag icon to reveal the name of the Coordinator that Flagged Member for follow-up.

Select the Members name to continue with Notifications workflow.