version:1.0



User Manual

Scale & Body Fat Analyzer Model # 21300





Contains FCC ID: REW-LSN8001

Manufactured for Zewa, Inc.

12960 Commerce Lakes Drive # 29 Fort Myers, FL 33913 USA www.zewa.com Toll Free Customer Service: 1-888-993-3592 warranty@zewa.com

- Thank you very much for selecting the Zewa Scale & Body Fat Analyzer 21300.
- Please do read the user manual carefully and thoroughly so as to ensure the safe usage of this product, and keep the manual well for further reference in case you have problems.

Table of Contents

Safety Information	
Safety and Usage Information	2
Your Scale and Its Environment	3
Efficient Use of Your Scale	3
Overview	
Device Components ·····	4
LCD Display	5
Initial Start-Up	
General Instructions	6
Insert the Batteries	6
Pair-up	7
Start Measuring	
Measurement Unit	8
Daily Measurement	8
Manage Your Health	
Data Transmission	9
Data Recognition	9
Troubleshooting	
Error Prompt	10
When Measuring	10-11
When Data Transmitting	12
Specifications	13
Maintenance	14
FCC Regulations	15
Appendix	
Table of Body Weight BMI	16
l able of Body Fat Level	16
Health Lips – About Body Fat	16
EMC Guidance	17-20

♥ Safety and Usage Information

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Please kindly find the meanings of the warning signs and symbols, which you may encounter in the label and user manual, as follows:

I	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "MANUFACTURER"
Bluetooth [*]	The Bluetooth Combination Mark	\sim	Symbol for "MANUFACTURE DATE"
★	Symbol for "TYPE BF APPLIED PARTS"		Symbol for "ENVIRONMENT
SN	Symbol for "SERIAL NUMBER"	X	PROTECTION – Waste electrical products should not be disposed of with
	Symbol for "DIRECT CURRENT"		household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "INDOOR USE ONLY"		

CAUTION

Zewa's Body Fat Analyzer 21300 offers you a seamless way to manage your health. Please be aware that this device is designed for adults' self-measuring and self-monitoring body fat level. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. This device should not be used by anyone who is acutely or chronically ill, suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.

21300 is equipped with data transmission function. It may emit electromagnetic energy so as to perform its intended function. Nearby portable and mobile RF communications equipment can affect the performance of 21300.

Portable and mobile RF communications equipment can affect the measuring accuracy of 21300.

Kindly note that the use of accessories, transducers or cables other than those specified, with the exception of transducers and cables sold by the manufacturer as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of 21300.

Be aware that misuse of electrical equipments can cause electric shock, burns, fire and other hazards.

Warning that the 21300 should not be used adjacent to or stacked with other equipment. Manufacturer will make available on request circuit diagrams, component parts list, etc. WARNING: No modifications of this equipment is allowed.

During using the Analyzer, the patient will contact with the surface and the electrode of the Analyzer.

INDICATIONS FOR USE

- The Zewa Body Fat Analyzer measures weight and uses bio-electrical impedance analysis (BIA) technology to estimate body fat in generally healthy adults 18 years of age or older.
- · It is intended for use in the home / domestic setting only.

CONTRAINDICATIONS

- 1. This device is contraindicated for any female subject who may be suspected of, or is pregnant. Besides provided inaccurate readings, the affects of this device on the fetus are unknown.
- 2. This device is contraindicated for any person who is connected to a wearable or implantable electronic device or instrument such as a pacemaker or defibrillator.

Your Scale and Its Environment

To ensure your safety as well as the service life of your scale, please avoid using the scale under the following circumstances:

- Concurrent use of this device and implantable medical electronic instruments, e.g. Cardiac Pacemaker
- Concurrent use of this device and wearable medical electronic instruments, e.g. electrocardiograph
- Concurrent use of this device and other medical electronic instruments for life support, e.g. mechanical heart
- · Slippery floor such as tile floor
- \cdot Jumping onto the platform immediately after bath or with wet hands
- · Near a cell phone or microwave oven

Avoid storage in the following locations:

- · Where there is water
- Where the device may be exposed to extreme temperatures, humidity, moisture, direct sunlight, dust, or salt air
- · Where there is risk of shock or drop
- · Where you store chemicals or full of corrosive gases
- \cdot Where in reach of the infants or children

Efficient Use of Your Scale

To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the scale on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Step onto the platform with bare feet. Stand still and keep full contact with the electrodes until the measurement is complete.
- Start measurement at least two hours after Getting up or Dinning.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dinning.
- Always start measurement in the same time slot and on the same scale located on the same flat, hard surface.

♥ Device Components





♥ List

- 1. Body Fat Analyzer 21300
- 2. Four AAA-size Batteries (1.5V each)
- 3. User Manual

♥ LCD Display



kg	Kilogram	8	Successful Bluetooth Connection
st Ib	Stone		Data transmitting / pending to transmit to wireless wellness system
lb	Pound		Low Battery

General Instructions

Zewa Body Fat Analyzer 21300 applies BIA (Bio-impedance Analysis) technology. A small amount of weak current flows through the human body so as to detect the bio-impedance and estimate body fat. The electrical current is small and may not be felt.

This BIA technology is cheap, safe, non-invasive, toxic-free and harmless. It also possesses the characteristics of simple operation and abundant information.

The current mentioned above is less than 1mA. However, please be aware that anyone with an wearable or implantable medical electronic instrument, such as a pacemaker, must avoid using this device.

The intended use of this device is for adult's indoor use only.

♥ Insert the Batteries

- · Open the battery door in the back of the scale.
- Insert the batteries (4 x 1.5V AAA) into the battery compartment according to the polarity indications marked inside the compartment.
 - * The digits "8888" will be shown on the LCD.



• Close the battery door and wait until the digits "0.0lb" are shown on the LCD.



CAUTION

• When the symbol
papears, the device will power off. Then you shall replace with a new set of batteries. Please replace all four batteries at the same time. Do NOT mix the old batteries with the new one.

Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage.
 Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries.

♥ Pair-Up

With the advanced Bluetooth 4.0 technology applied, the mobile or portable equipments, which are equipped with Bluetooth function in line with BLE Technical Specifications established by global organization Bluetooth SIG, are capable to receive your personal health data.

- The Zewa App is now available in App Store
 You may search and install the Zewa app in your smartphone.
- · Turn on Bluetooth and the Zewa App. Make sure both are ON when pair-up is proceeding.
- Press and hold "UNIT" button in the back of the scale to start pair-up.



Bluetooth Module No. : nRF8001						
Frequency Range 2.402 - 2.480 GHz Supply Voltage 3.3 V						
Output Power Range	0 dBm					

Select Measurement Unit

With batteries correctly installed, press "UNIT" button in the back of the scale to select measurement unit. The default measurement unit is "pound". You may press "UNIT" button to choose among pound, kilogram and stone.



Daily Measurement

• With original SENSE ON patent technology, 21300 will automatically switch on as you step on the platform barefooted.



 Stand still and keep full contact with the electrodes until the LCD stops displaying moving "o".



 When your scale is successfully paired with your smartphone, please open the Zewa App and turn on your Bluetooth. Keep smartphone and scale in transmission distance. Then the measurement data will automatically transferred to the app.

(Please refer to <u>Data Transmission</u> for more details.)

♥ Data Transmission

Once 21300 is successfully paired-up with your smartphone, please open the Zewa App and turn on your Bluetooth. Keep smartphone and scale in transmission distance. Step on the scale and take a measurement. Once the measurement is finished the data is automatically transferred to the app.



- The symbol 🗐 will disappear after successful data transmission, and you may check your personal health data stored in your smartphone.
- If the data transmission fails, the symbol is will remain. The pending measurement data will be temporarily kept in the scale and transmitted to your smartphone when next measurement is complete.



_ 🤥 CAUTION _

- Interference may occur in the vicinity of equipment marked with the following symbol ((2)). And the Analyzer may interfere vicinity electrical equipment.
- To enable the data transmission function, this product should be paired to a Bluetooth 4.0 end.

How to mitigate possible interference?

- The range between the Analyzer and the Bluetooth end should be reasonably close, from 1 meter to 10 meters. Please ensure no obstacles between the Analyzer and the Bluetooth end so as to obtain quality connection.
- 2. To avoid interference, other electronic devices (particularly those with Bluetooth transmission / Transmitter) should be kept at least 1 meter away from the Analyzer.

Data Recognition

- After the weight data is transmitted successfully, you can check the measurement result in your smartphone.
- According to the results, the weight will match the claimed user with closest weight stored in the Zewa App automatically. If there are two or more users sharing very similay body characteristics, or if there is a big deviation between the connect user and the result, the App will require confirmation to assign the result, then you may select a user manual.

Troubleshooting

• Error Prompt

Error	Description	Solution
·····	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.
EI	Failure of pairing up.	Please check below items: -Bluetooth is ON. -App Collector is ON. -Both devices are within the transmission distance of Bluetooth.

♥ When Measuring ...

Problem	Root Cause	Solution
Abnormal measuring results: - Too high; OR - Too low; OR - Huge difference between two recent measurement.	Incorrect posture	Please step on the platform barefooted and stand still.
	The device is located on the soft ground such as a carpet OR on a rugged surface.	Please place the device on a flat, hard surface.
	Cold body that may results in bad blood circulation.	Warm up your hands and feet to resume blood circulation and then measure again.

Problem	Root Cause	Solution
Abnormal measuring results: - Too high; OR	Cold Electrodes.	Place the device in a warm room for a while and then measure again.
- Huge difference between two recent measurement.	Either your hands or your feet are too dry.	Wipe your feet with a damp cloth, keeping them slightly damp when starting measurement.
No display on LCD when the device powers on.	Batteries not yet installed.	Install the batteries. (Please refer to <u>Insert the</u> <u>Batteries</u>)
	Worn batteries.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.
CANNOT proceed to analyze body fat.	Step onto the platform wearing socks or shoes.	Please keep barefooted during the measurement, and keep full contact with the electrodes as well.
The device powers off automatically.	Low battery.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.

♥ When Data Transmitting ...

Problem	roblem Root Cause Solution			
Data transmission failed.	Bluetooth is OFF.	Turn ON the Bluetooth via "Setting >> General >> Bluetooth".		
	App is OFF.	Press the icon to turn ON your app.		
	Out of range of Bluetooth transmission.	Place your smartphone closer to the scale.		

♥ Specifications

Product Name	Body Fat Analyzer (BF-1255-B)				
Dimension	Scale: 310x310x31mm				
Net Weight	Approxim ately 1.63kg (Excluding the dry cells)				
Display	Digital LCD with White Backlight				
Measurement Unit	Kilogram / Stone/ Pound				
Measurement Range	5kg to 180kg / 0st: 11lb to 28st: 5lb / 11lb to 397lb				
Division	0.1kg / 0.2lb				
Accuracy	5-50kg: ±0.3kg; 50-100kg: ±0.4kg; 100-150kg: ±0.5kg; 150-180kg: ±0.7kg				
Working Environment	Temperature: 0 ℃ to 40 ℃ Humidity: ≤90% RH				
Storage Environment	Temperature: -20 C to 60 C Humidity: 10%RH to 93% RH				
Power Source	6V (Four AAA-size Batteries)				
Turn on Method	SENSE ON technology				
Auto-OFF	The scale will turn off after about 10 seconds if there is no operation				
Accessories	1. Four AAA-size batteries 2. User Manual				
Mode of Operation	Continuous Operation				

About the Accuracy of This Product

• This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer. Please refer to the above table for the descriptions on accuracy.

• This product is specially designed for body fat analysis as well as weight measurement. It should NOT be used by anyone during the process of transaction for verification of goods' weight.

Maintenance

♥ Maintenance

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems, please contact Transtek. (*Please refer to <u>Warranty</u> for contact information*)
- Do not dispose of batteries in fire.Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period.

♥ FCC Regulations

FCC User Guide Information

Radio Frequency Interface Requirements - FCC

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving antenna;
- · Increase the separation between the equipment and receiver;
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected;
- · Consult the dealer or an experienced radio / TV technician for help.

Radio Transmitters (Part 15)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC RF Exposure Guidelines

Safety Information

Reducing RF Exposure - Use Properly

Only operate the device in accordance with the instructions supplied. This device complies with FCC radiation exposure limits set forth for an uncontrolled environment. In order to avoid the possiblity of exceeding the FCC radio frequency exposure limits, human proximity to the antenna shall not be less than 20 cm (8 inches) during normal operation.

Appendix

♥ Table of Body Weight BMI

BMI=weight (kg) / [height (m)] ² BMI=weight (lb) / [height (in)] ² ×703

Rating	Range	Colour
Underweight	BMI<18.5	Blue
Normal	18.5≤BMI≤24.9	Green
Overweight	25≤BMI≤29.9	Orange
Obesity	BMI≥30	Red

♥ Table of Body Fat Level (Unit: %)

a) The body fat percentage (%): 5%-60%/0.1%

Standard for Men

Standard for Women

Rating			Age			Rating			Age	50.50	
	20-29	30-39	40-49	50-59	60+		20-29	30-39	40-49	50-59	-00+
low	<13	<14	<16	<17	<18	low	<19	<20	<21	~22	~23
Normal	14-20	15-21	17-23	18-24	19-25	Normal	20-28	21-29	22-30	23-31	24-52
Moderately High	21-23	22-24	24-26	25-27	26-28	Moderately High	29-31	30-32	31-33	32-33	33-35
High	>23	>24	>26	>27	>28	High	>31	>32	>33	204	- 55

Source: University of Illinois Department of Food Science and Human Nutrition. Body Fat Percentage Calculator.

www.ag.uiuc.edu/~food-lab/ai/bfc.html

♥ Health Tips - About Body Fat

Fat is essential for human body. It can not only store energy and protect viscera, but also regulate body temperature and maintain normal physiological function of human body. However, too much body fat is harmful to human body. It is always accompanied by Fatty Liver, diabetes, coronary heart disease, etc.

Therefore self-measuring and self-monitoring body fat level are beneficial to your health. Since we can't judge body fat level simply by our weight, this body fat analyzer 21300, with BIA (Bio-impedance Analysis) technology applied, is an accurate device that offers a quick and comfortable way to obtain your body fat level.

• EMC Guidance

Table 1 – Guidance and MANUFACTURER'S declaration – ELECTROMAGNETIC EMISSIONS for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions					
The device is intended The customer or the u environment.	The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.				
Emissions test	Compliance	Electromagnetic environment-guidance			
RF emissions CISPR 11	Group 2	The device must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.			
RF emissions CISPR 11	Class B				
Harmonic emissions IEC 61000-3-2	Not applicable				
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable				

Table 2 – Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity						
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.						
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance			
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.			
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.			
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.			
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U ₇ (95 % dip in U ₇) for 0,5 cycle 40 % U ₇ (60 % dip in U ₇) for 5 cycles 70 % U ₇ (30 % dip in U ₇) for 25 cycles <5 % U ₇ (95 % dip in U ₇) for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uniterruptible power supply or a battery.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.			
NOTE U_T is the a.c. mains voltage prior to application of the test level.						

Table 4 – Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING



Appendix

Table 6 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Beted maximum autority and	Separation distance according to frequency of transmitter m			
of transmitter	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz	
w	$d = [\frac{3,5}{V_1}]\sqrt{P}$	$d =$ 1.167 \sqrt{P}	$d = 2.333 \sqrt{P}$	
0,01	Not applicable	0.117	0.233	
0,1	Not applicable	0.369	0.738	
1	Not applicable	1.167	2.333	
10	Not applicable	3.690	7.378	
100	Not applicable	11.67	23.33	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.