

Pairing the Acon Glucose Meter with the Life365 App

 Follow Manufacturer Instructions for Use.
For additional assistance, contact your Care Provider.

Version Update August 2020

- 1 Select the **Hamburger Icon (☰)**, located in the upper left of the Life365 app screen. Select **My Devices**.
 - 2 On the My Devices screen, select the button **Add New**. The smart device will search for your Bluetooth device.
 - 3 Turn **ON** the ACON OnCall Express Mobile Glucose Meter by holding the **OK** button for 3-4 seconds. *Complete initial setup of the Glucose meter (Date, Time etc.).*
 - 4 On the glucose meter screen, the following text will display: **“no PAIR bt”**

Press the up-arrow once until the glucose meter screen displays **“YES PAIR bt”**. Select the **OK** button.

A small green light on the right side of the device will flash while the device attempts to connect.
 - 5 Once found in the app, select the Bluetooth device **name**.
 - 6 Select the button **Add to My Devices**.
 - 7 A prompt may display requesting a code to complete device pairing. If a prompt does not display, complete the following:
 - With an index finger, swipe downwards from the top of your smart device screen to reveal Notifications.
 - Select the notification **Pairing Request**.
 - Enter the **6-digit code** displayed on the Glucose meter screen.
- Check the box **“Allow access to your contacts and call history”**. Select **OK** once completed.

Device pairing is now complete.

For video instructions, and other helpful resources, visit
<https://www.life365.health/en/support/app>.