



## QUICK PAIRING GUIDE

# Pairing the TaiDoc / ForaCare Thermometer with MyLife365 Smart Devices

 Follow Manufacturer Instructions for Use.

*Still have questions or need assistance? Contact  
Support at [support@life365.health](mailto:support@life365.health)*

# Pairing the TaiDoc / ForaCare Thermometer with MyLife365 Smart Devices

1

After assigning a smart device to a member, select the Reminders application. You will receive two prompts, select “**ALLOW**” for both prompts.

**This step must be completed before pairing Bluetooth Devices.**



Allow **Reminders** to make and manage phone calls?

1 of 2

DENY ALLOW



Allow **Reminders** to access photos, media, and files on your device?

2 of 2



DENY ALLOW

2

Select the Bluetooth Icon\*.

Enter the password “**3651**” and select “**OK**”.



*\*the Bluetooth App Icon location may vary by device*



3

On the “**My Devices**” page, select “**Add New**” (A).

The smart device will search for the Bluetooth device.



# Pairing the TaiDoc / ForaCare Thermometer with MyLife365 Smart Devices

4



The TaiDoc/Fora Thermometer requires a measurement reading to complete initial pairing:

- 4.1 Insert an ear protection cover on the metal tip of the thermometer.
- 4.2 Press the power button to turn ON the thermometer.
- 4.3 Take a measurement reading with the thermometer.

Allow a few moments for the Bluetooth pairing application to find the thermometer.

5


Once found, select the Bluetooth device **name** (B).

Select “Add to My Devices” (C).

*This is an example. Device name displayed will reflect the Bluetooth device the user is attempting to pair.*

6

Device pairing is now complete.




# MANUAL PAIRING:

## Difficulty pairing this device?

- 1 Try manually pairing the device by selecting “**Add New**” on the “**My Devices**” page (A), then selecting “**Manually**” (B).
- 2 Search for the device by Manufacturer by selecting the magnifying glass (C).
- 3 Select the device, then select “**Add to My Devices**”.

Device pairing is now complete.

Take a test reading with the device.  
Hub will announce the reading.



### NOTE:

The Reminders Application must be opened BEFORE opening the Bluetooth Application. If the Bluetooth Application is opened before the Reminders Application, a “**Choose Server Location**” error screen will display.

If this error screen is shown:

1. Exit the Bluetooth Application
2. Enter the Reminders Application
3. Accept both prompts
4. Re-enter the Bluetooth Application
5. Select the Refresh symbol (D).

**You can now complete device pairing.**

