

# Product Data Sheet



## ProFoam Positioning Aid

Ref.:	Product name:	Dimensions:	Weight:	Packaging dimensions:
1150-B	ProFoam Plate 25x18x5	250 x 180 x 50 mm	120 g	280 x 470 mm
1151-B	ProFoam Wedge 18x25x12	180 x 250 x 120 mm	160 g	350 x 500 mm
1152-B	ProFoam Wedge 26x25x7	260 x 250 x 70 mm	165 g	350 x 500 mm
1153-B	ProFoam Wedge 18x11x5	180 x 110 x 50 mm	40 g	280 x 470 mm
1154-B	ProFoam Wedge 36x25x12	360 x 250 x 120 mm	275 g	350 x 500 mm
1155-B	ProFoam Wedge 60x25x12	600 x 250 x 120 mm	440 g	500 x 800 mm
1156-B	ProFoam Plate 32x32x10 DS	320 x 320 x 100 mm	785 g	500 x 800 mm
1157-B	ProFoam Plate 32x32x5 DS	320 x 320 x 50 mm	445 g	400 x 600 mm
1158-B	ProFoam Plate 25x18x5 DS	250 x 180 x 50 mm	195 g	280 x 470 mm
1159-B	ProFoam Wedge 18x25x12 DS	180 x 250 x 120 mm	250 g	350 x 500 mm
1160-B	ProFoam Wedge 26x25x7 DS	260 x 250 x 70 mm	225 g	350 x 500 mm
1161-B	ProFoam Wedge 18x11x5 DS	180 x 110 x 50 mm	60 g	280 x 470 mm
1162-B	ProFoam Wedge 36x25x12 DS	360 x 250 x 120 mm	470 g	350 x 500 mm
1163-B	ProFoam Wedge 60x25x12 DS	600 x 250 x 120 mm	785 g	500 x 800 mm
1181-A	ProFoam Plate DIM	variabel	variabel	variabel
1182-A	ProFoam Wedge DIM	variabel	variabel	variabel

<b>Intended use</b>	Short-term patient positioning and stabilization device during invasive and non-invasive procedures to be used on intact patient's skin. Non-Sterile and not to be used in the sterile field.	
<b>Directives &amp; classification</b>	CE according to MDD – Directive 93/42/ECC and MDR – Regulation 2017/745 Class 1 product	
<b>Applied quality assurance</b>	SN EN ISO 13485:2016	
<b>Instructions for use</b>	Enclosed to the product	
<b>Restrictions of use</b>	Do not use the product for purposes other than those described.	
<b>General safety notes</b>	Non sterile. Do not sterilize. Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).	
<b>MR safety note</b>	MR Safe	
<b>Materials</b>	<b>Material</b> TPU PU (Polyether based) LDPE	<b>Component</b> foils foam packaging
<b>Marking</b>	CE	
<b>Storage</b>	Store in dry place / room temperature	

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.