

Product Data Sheet



ProBelt

Ref.:	Product name:	Dimensions:	Weight:	Packaging dimensions:
1170-A	ProBelt 200 R2	200 x 650 mm each side	240 g	400 x 270 x 20 mm
1171-A	ProBelt 200 R3	200 x 650 mm each side	240 g	400 x 270 x 20 mm
1172-A	ProBelt 200 T5	200 x 650 mm each side	240 g	400 x 270 x 20 mm
1173-A	ProBelt Thorax 200 PFC	Cushion: 500 x 300 x 50 mm Belt buttons: 400 x 200 mm Belt buttonholes: 600 x 200 mm	300 g	400 x 600 x 50 mm
1174-A	ProBelt Thorax 200 PFW	Wedge: 300 x 500 x 150 mm Belt buttons: 400 x 200 mm Belt buttonholes: 600 x 200 mm	700 g	500 x 800 x 120 mm

Intended use	Short-term patient positioning and stabilization device during invasive and non-invasive procedures to be used on intact patient's skin. Non-Sterile and not to be used in the sterile field.	
Directives & classification	CE according to MDD – Directive 93/42/ECC and MDR – Regulation 2017/745 Class 1 product	
Applied quality assurance	SN EN ISO 13485:2016	
Instructions for use	Enclosed to the product	
Restrictions of use	Do not use the product for purposes other than those described.	
General safety notes	Non sterile. Do not sterilize. Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).	
MR safety note	MR Safe	
Materials	Material TPU TPU/ Nylon PU/ PVC PVC LDPE	Component foils textile foil buttons profile packaging
Marking	CE	
Storage	Store in dry place / room temperature	

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.