

Pearltec Positioning Aids

PearlFit Vac

PearlFit Vac positioning aids become stable by withdrawing air and allow a reproducible patient positioning.



PearlFit Cushion Vac

30 x 20 x 4 cm Vac	Art. 1125
50 x 20 x 4 cm Vac	Art. 1131
50 x 30 x 7 cm Vac	Art. 1126
50 x 30 x 10 cm Vac	Art. 1127
50 x 50 x 5 cm Vac	Art. 1143
90 x 60 x 5 cm Vac	Art. 1144
130 x 60 x 10 cm Vac	Art. 1129
200x 75 x 7 cm Vac	Art. 1141

Depending on the size, PearlFit Cushions Vac are perfectly suitable for a stable patient positioning inside and outside the field of view without any pressure points.



PearlFit Wedge Vac

26 x 25 x 7 cm Vac	Art. 1128
30 x 50 x 15 cm Vac	Art. 1139
56 x 50 x 25 cm Vac	Art. 1133

The PearlFit Wedge Vac helps to create a specific angle for certain body parts.



PearlFit Head/Neck Vac

Art. 1130

The PearlFit Head/Neck Vac is used to support the patient's head when lying in supine position.



Hand Vacuum-Pump

Art. buw-5355

The Hand Vacuum-Pump is MR-compatible and can be used to withdraw air from vacuum cushions.

