


Pearltec Positioning Aids

PearlFit


PearlFit positioning aids can be used for stable but yet comfortable patient positioning inside and outside the field of view during all radiological procedures.



PearlFit Cushion

22 x 8 x 1 cm		Art. 1136
27 x 6 x 1 cm		Art. 1145
30 x 11 x 2 cm		Art. 1137
30 x 20 x 4 cm	TIP	Art. 1114
50 x 20 x 4 cm		Art. 1124
50 x 30 x 7 cm	TIP	Art. 1115
50 x 30 x 10 cm	TIP	Art. 1116
50 x 50 x 5 cm	TIP	Art. 1142
130 x 60 x 10 cm		Art. 1120
200x 75 x 7 cm		Art. 1140

Depending on the size, PearlFit Cushions are perfectly suitable for a stable patient positioning inside and outside the field of view without any pressure points.



PearlFit Wedge

26 x 25 x 7 cm		Art. 1117
30 x 50 x 15 cm		Art. 1138
56 x 50 x 25 cm	TIP	Art. 1132

The PearlFit Wedge helps to create a specific angle for certain body parts.



PearlFit Roll

Ø = 18, L = 40 cm		Art. 1112
Ø = 18, L = 50 cm	TIP	Art. 1113

The PearlFit Roll can be used during all radiological procedures and is primarily used in the hollow of the knee when patient is lying in supine position.



PearlFit Flex

6 chambers, 90 x 50 x 5 cm

Art. 1118

3 chambers, 90 x 50 x 7 cm

Art. 1119

The PearlFit Flex products can be used as either cushion, wedge or roll making them universally applicable for a broad range of examinations.



PearlFit Head/Neck

TIP

Art. 1121

The PearlFit Head/Neck is used to support the patient's head when lying in supine position.



PearlFit Forehead

Art. 5005

The PearlFit Forehead consists of two components and is used to support the patient's head in prone position.



PearlFit Inlay Head

Art. 1123

PearlFit Inlay is form-fitted to a broad range of coils, cradles and holders and reduces pressure points at the back of the head.