


Pearltec Positioning Aids

ProFoam – soft


The ProFoam positioning aids are made of a soft foam material and are welded into a durable polyurethane foil. They are available in the most common shapes and sizes and are suitable for a wide range of applications, such as in the field of radiology.



ProFoam Plate **NEW**

25 x 18 x 5 cm


Art. 1150



ProFoam Wedge 45 **NEW**

18 x 25 x 12 cm
18 x 11 x 5 cm
36 x 25 x 12 cm
60 x 25 x 12 cm

Art. 1151
Art. 1153
Art. 1154
Art. 1155



ProFoam Wedge 15 **NEW**

26 x 25 x 7 cm


Art. 1152

Other shapes and sizes available upon request.

Pearltec Positioning Aids

ProFoam – stable/hard


The ProFoam positioning aids DS are made of a dimensionally stable foam material and are welded into a durable polyurethane foil. They are available in the most common shapes and sizes and are suitable for a wide range of applications, such as in the field of radiology.



ProFoam Plate DS **NEW**

25 x 18 x 5 cm


Art. 1158



ProFoam Wedge 45 DS **NEW**

18 x 25 x 12 cm
18 x 11 x 5 cm
36 x 25 x 12 cm
60 x 25 x 12 cm

Art. 1159
Art. 1161
Art. 1162
Art. 1163



ProFoam Wedge 15 DS **NEW**

26 x 25 x 7 cm

Art. 1160

Other shapes and sizes available upon request.