

iFIX Fixation System

The iFIX Patient Fixation System is a hygienic alternative to belts, tapes and hook and loop. It combines micro-hook anchor products with a medical single-use fleece.



iFIX Fleece

5 (Roll: 5 cm x 30 m)
15 (Roll: 15 cm x 200 m)
50 (Roll: 50 cm x 30 m)

Art. iSYS-31618
Art. iSYS-16304
Art. iSYS-17674

iFIX Fleece is available in three different sizes and comes on a roll in a dispenser box. By utilizing the included cutter it can be cut to the length needed.



iFIX Head Mask

XS (10 pcs)
S (10 pcs)
M (10 pcs)
L (10 pcs)
XL (10 pcs)
L elongated (10 pcs)

Art. iSYS-26102
Art. iSYS-26105
Art. iSYS-26104
Art. iSYS-26103
Art. iSYS-29341
Art. iSYS-35080

iFIX Head Mask is used to immobilize the patient's head.



iFIX Slot Adapter

8/30 (4 pcs)
3/30 (4 pcs)

Art. iSYS-16192L003
Art. iSYS-17622L003

iFIX Slot Adapter is coated with micro-hook and can be slid into the rails on the patient table.



iFIX Patch

15 (8 pcs)
25 (6 pcs)
50 (4 pcs)
75 (4 pcs)
15 x 6 (3 pcs)

Art. iSYS-16180L001
Art. iSYS-16180L002
Art. iSYS-16180L005
Art. iSYS-16180L007
Art. iSYS-32982

iFIX Patch is a self-adhesive micro-hook that can be applied on the patient table as well as on PearlFit (Vac) positioning aids cushions.



iFIX Patch Adapter

50 (4 pcs)

Art. iSYS-16179L001

iFIX Patch Adapter is coated with micro-hook and can be applied on existing hooks of holders, cradles and patient tables.

