



The wait is over.

Introducing Atmos—the first innovation in air-charged urodynamics catheters in twenty years. Designed for accurate results, efficient usability, and maximum comfort, we invite you to see why physicians are choosing Atmos.

A better choice for physicians and their patients



Atmos¹ Air-Charged Urodynamics Catheters are designed to match today's leading catheter while bringing added enhancements to improve the diagnostic experience for patients and practitioners.

- Flexible design for a relaxed fit
- The smallest diameter on the market
- Compatible with your existing equipment
- Consistent and reliable results

Proven accuracy and performance

Tested with 151 catheters at 1,991 pressure points, Atmos delivered near perfect statistical agreement—even in equipment calibrated to a T-DOC® Air-Charged™ Urodynamics Catheter.²

STATISTICAL APPROACH	STATISTIC	95% CONFIDENCE INTERVAL	P
Lin's Concordance Correlation Coefficient P _c	.999	.999, >.999	<.001
Pearson's r	1.00		<.001
Bland & Altman's Limits of Agreement	-0.897	-2.884, 1.090	



Everything you need in one place.

We provide an assortment of quality essential urodynamic consumable items to simplify procurement creating efficiency for your practice.

Consumables

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	BOX QTY
1A-AC	Single Sensor Abdominal Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1B-AC	Single Sensor Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1CB-AC	Single Sensor Coudé Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
2BU-AC	Dual Sensor Bladder/Urethral Air-Charged Urodynamics Catheter	6 FR, 25" long	10

Pump Tubing and Extension Lines

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	BOX QTY
UD-A36205	Pump Tubing, 3,2x6,4 mm	Laborie Legacy Equipment	25
UD-A58245	Pump Tubing, 5x8	Medical Measurement System Equipment	25
MMT-1173	Urodynamic Set for Pump	MEDKONSULT Medical Technology Equipment	25
UD-A46215	Pump Tubing, 4x6 mm	Mediwatch Equipment	100
UD-AELN	Transparent Extension Line, 1500 mm	-	50

Electrodes

PRODUCT NUMBER	DESCRIPTION	BOX QTY
CND1720003 - 10	Cleartrode Repositionable ECG Electrodes Snap Patch (10 Sets)	30
CND1720003 - 100	Cleartrode Repositionable ECG Electrodes Snap Patch (100 Sets)	300
RA36-03	Three 36"/91.4cm Snap Leadwires (Black, Green, White)	3
RA36-03II	Three 36"/91.4cm Snap Leadwires (Black, Red, White)	3
CND1741003 - 10	Neotrode II Pediatric/Neonatal ECG Electrode, Tape with Wire Attached, 3/Pouch (10 Sets)	30
CND1741003 - 100	Neotrode II Pediatric/Neonatal ECG Electrode, Tape with Wire Attached, 3/Pouch (100 Sets)	300

Choose Atmos.

Finally, you have access to premium products at a responsible price and can choose what's best for your practice.

Visit geminimedtech.com to request a sample or place an order.

GEMINI

Gemini Medical Technologies LLC
123 NW 13th Street, Suite 212
Boca Raton, FL 33432

Phone: +1.833.943.6464
Email: info@geminimedtech.com
geminimedtech.com

* GEMINI and ATMOS are trademarks of Gemini Medical Technologies LLC. All Rights Reserved. Product names, brands and other trademarks featured or referred to within the GEMINI website and marketing documents are the property of their respective trademark owners

* T-DOC® is a registered trademark of Laborie Medical Technologies Corp

