

# YOU PRESCRIBE BIOFINITY® TORIC SO MUCH, WE MADE IT TWICE.


The same optical design of the #1, most prescribed toric lens on the market<sup>1</sup> is available in MyDay® toric. The silicone hydrogel lens brings consistent centration, stable fit and excellent visual acuity that you and your patients can rely on.

T O R I C **2** W A Y S



## UNMATCHED STABILITY<sup>2</sup>

Optimized Toric Lens Geometry™




creates a multifaceted toric lens that's designed to provide predictable, consistent visual acuity, lens stability, fit and comfort as well as better performance and simple fitting through:

- Uniform horizontal iso-thickness
- Smooth, continuous ballast design
- Large toric optic zone

## UNSURPASSED COMFORT<sup>3</sup>

Aquaform® Technology



provides natural moisture and breathability:

- Long silicone chains and compatible hydrogen bonds lock water molecules within the lens
- Less silicone is needed to optimize oxygen transmissibility
- Reduced silicone content results in an optimum modulus, making the lenses soft and flexible

# IT'S A MONTHLY. IT'S A DAILY. IT'S TORIC 2 WAYS.



**MyDay® toric** lenses ensure 100% corneal oxygen consumption, transmitting high levels of oxygen across the entire lens surface to help maintain ocular health<sup>4</sup>. They also filter 86% of UVA and 97% of UVB rays<sup>5</sup>. Plus, MyDay® toric has the most prescription options of any 1-day lens.<sup>6</sup>

<b>Material/H2O content</b>	stenfilcon A/54%
<b>Replacement schedule</b>	Daily
<b>Oxygen transmissibility</b>	80 Dk/t (at -3.00D)
<b>Revenue carton size</b>	90-pack
<b>Design</b>	Optimized Toric Lens Geometry™
<b>Base curve</b>	8.6 mm
<b>Diameter</b>	14.5 mm
<b>UV Blocker</b>	Yes

Sphere power	Cylinder	Axis
Plano to -6.00D (0.25D steps)	-0.75, -1.25, -1.75	Full circle in 10° steps
Plano to -6.00D (0.25D steps) -6.50D to -10.00D (0.50D steps) +0.50D to +6.00D (0.50D steps)	-2.25	10°, 20°, 90°, 160°, 170°, 180°
-6.50D to -10.00D (0.50D steps) +0.25D to +6.00D (0.25D steps)	-0.75, -1.25, -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°

**Biofinity® toric** and **Biofinity® XR toric** have more prescription options than all other monthly replacement SiHy toric brands combined<sup>7</sup>. Sphere powers expanded to ±20.00D bring the total number of unique prescription options to over 37,000.

<b>Material/H2O content</b>	comfilcon A/48%
<b>Replacement schedule</b>	Monthly
<b>Oxygen transmissibility</b>	160 Dk/t (at -3.00D)
<b>Extended wear</b>	Yes
<b>Revenue carton size</b>	6-pack blisters
<b>Design</b>	Asphere
<b>Base curve</b>	8.6 mm
<b>Diameter</b>	14 mm

Sphere power	Cylinder	Axis
+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75, -1.25, -1.75, -2.25	10° to 180° (in 10° steps)
+8.50D to +20.00D and -10.50D to -20.00D (0.50D steps)	-0.75, -1.25, -1.75, -2.25	5° to 180° (in 5° steps)
+20.00D to -20.00D (0.50D steps after +/-6.00D)	-2.75, -3.25, -3.75, -4.25, -4.75, -5.25, -5.75	5° to 180° (in 5° steps)

With our OptiExpert™ fitting tool app, speed up finding the right toric contact lens fit for your patients.

**Contact your CooperVision sales representative or Authorized Distributor to get started.**