


SPECIFICATIONS



Transonic® Clinical Flowprobes

Models:

AU, FD, FMC, FME, FMV, MB/MR, MB-S/MR-S, MU, P

transonic
THE MEASURE OF BETTER RESULTS.

USA/Canada



Transonic Systems Inc.
34 Dutch Mill Road
Ithaca, NY 14850
U.S.A.
Tel: +1 607-257-5300
Fax: +1 607-257-7256
support@transonic.com



Transonic Europe B.V.
Business Park Stein 205
6181 MB Elsloo
The Netherlands
Tel: +31 43-407-7200
Fax: +31 43-407-7201
europe@transonic.com

Asia/Pacific

Transonic Asia Inc.
6F-3 No 5 Hangsiang Rd.
Dayuan, Taoyuan County
33747 Taiwan, R.O.C.
Tel: +886 3399-5806
Fax: +886 3399-5805
support@transonicasia.com

Japan

Nipro-Transonic Japan Inc.
7th Floor, Maruha Building
11-1 Matsuba-cho
Tokorozawa City, Saitama
359-0044 Japan
Tel: +81 04-2946-8541
Fax: +81 04-2946-8542
japan@transonic.com



Clinical Flowprobe Specifications

Probe Size	Catalog #	Max Flow	Accuracy	Max Zero Offset (ml/min)		Ultrasound Frequency
mm		L/min	% of Reading	HQC/D/N	HQE/F/G	MHz
FMC-Series: Coronary immediate-use steam sterilizable Probes (long handle, J-reflector)						
Available with HQC, HQD or HQF prefix						
1.5mm	HQ_1.5FMC	0.5	+/- 10	+/- 2	+/- 1	4.8
2mm	HQ_2FMC	0.5	+/- 10	+/- 2	+/- 1	4.8
3mm	HQ_3FMC	1	+/- 10	+/- 3	+/- 1.5	3.6
4mm	HQ_4FMC	2	+/- 10	+/- 6	+/- 3	2.4
FMV-Series: Vascular immediate-use steam sterilizable Probes (short handle, J-reflector)						
Available with HQC, HQD or HQF prefix						
1.5mm	HQ_1.5FMV	0.5	+/- 10	+/- 2	+/- 1	4.8
2mm	HQ_2FMV	0.5	+/- 10	+/- 2	+/- 1	4.8
3mm	HQ_3FMV	1	+/- 10	+/- 3	+/- 1.5	3.6
4mm	HQ_4FMV	2	+/- 10	+/- 6	+/- 3	2.4
6mm	HQ_6FMV	5	+/- 10	+/- 15	+/- 7.5	1.8
8mm	HQ_8FMV	10	+/- 10	+/- 30	+/- 15	1.2
10mm	HQ_10FMV	10	+/- 10	+/- 30	+/- 15	1.2
12mm	HQ_12FMV	20	+/- 10	+/- 60	+/- 30	0.9
14mm	HQ_14FMV	20	+/- 10	+/- 60	+/- 30	0.9
FME-Series: Carotid Endarterectomy immediate-use steam sterilizable Probes (short handle, L-reflector)						
Available with HQC*, HQD, HQN* or HQF prefix						
1.5mm	HQ_1.5FME	0.5	+/- 10	+/- 2	+/- 1	4.8
2mm	HQ_2FME	0.5	+/- 10	+/- 2	+/- 1	4.8
3mm	HQ_3FME	1	+/- 10	+/- 3	+/- 1.5	3.6
4mm	HQ_4FME	2	+/- 10	+/- 6	+/- 3	2.4
6mm	HQ_6FME	5	+/- 10	+/- 15	+/- 7.5	1.8
8mm	HQ_8FME	10	+/- 10	+/- 30	+/- 15	1.2
10mm	HQ_10FME	10	+/- 10	+/- 30	+/- 15	1.2
MB-Series: Intracranial Charbel Micro-Flowprobes (long handle, single use)						
Available with HQN or HQG prefix						
1.5mm	HQ_1.5MB	0.2	+/- 10	+/- 3	+/- 1.5	4.8
2mm	HQ_2MB	0.5	+/- 10	+/- 5	+/- 2.5	3.6
3mm	HQ_3MB	1	+/- 10	+/- 10	+/- 5	2.4
MB-S-Series: Extracranial Micro-Flowprobes (short handle, single use)						
Available with HQN or HQG prefix						
3mm	HQ_3MB-S	1	+/- 10	+/- 3	+/- 1.5	3.6
4mm	HQ_4MB-S	2	+/- 10	+/- 6	+/- 3	2.4
6mm	HQ_6MB-S	5	+/- 10	+/- 15	+/- 7.5	1.8
MR-Series: Intracranial Charbel Micro-Flowprobes (long handle, multi-use)						
Available with HQN or HQE prefix**						
1.5mm	HQ_1.5MR	0.2	+/- 10	+/- 3	+/- 1.5	4.8
2mm	HQ_2MR	0.5	+/- 10	+/- 5	+/- 2.5	3.6
3mm	HQ_3MR	1	+/- 10	+/- 10	+/- 5	2.4

Clinical Flowprobe Specifications

Probe Size	Catalog #	Max Flow	Accuracy	Max Zero Offset (ml/min)		Ultrasound Frequency
mm		L/min	% of Reading	HQC/D/N	HQE/F/G	MHz
MR-S-Series: Extracranial Micro-Flowprobes (short handle, multi-use) Available with HQN or HQE prefix**						
3mm	HQ_3MR-S	1	+/- 10	+/- 3	+/- 1.5	3.6
4mm	HQ_4MR-S	2	+/- 10	+/- 6	+/- 3	2.4
6mm	HQ_6MR-S	5	+/- 10	+/- 15	+/- 7.5	1.8
P-Series: Port Access Probes (long handle, J-reflector) Available with HQF prefix***						
2mm	HQ_2P	0.5	+/- 10	N/A	+/- 1	4.8
3mm	HQ_3P	1	+/- 10	N/A	+/- 1.5	3.6
4mm	HQ_4P	2	+/- 10	N/A	+/- 3	2.4
FD- Series: Port Access Probes with Couplant delivery set (long handle, J-reflector) Available with HQF prefix						
1.5mm	HQ_1.5FD	0.5	+/- 10	N/A	+/- 0.8	4.8
2mm	HQ_2FD	0.5	+/- 10	N/A	+/- 1	4.8
3mm	HQ_3FD	1	+/- 10	N/A	+/- 1.5	3.6
4mm	HQ_4FD	2	+/- 10	N/A	+/- 3	2.4
MU-Series: Microvascular Probes (short handle, L-reflector) Available with HQF prefix						
0.7mm	HQ_0.7MU	50 mL/min	+/- 10	N/A	+/- 0.4	9.6
1mm	HQ_1MU	100 mL/min	+/- 10	N/A	+/- 0.8	7.2
1.5mm	HQ_1.5MU	200 mL/min	+/- 10	N/A	+/- 1.5	4.8
2mm	HQ_2MU	500 mL/min	+/- 10	N/A	+/- 2.5	3.6
3mm	HQ_3MU	1000 mL/min	+/- 10	N/A	+/- 5	2.4
AU-Series: Confidence Flowprobes Size 8-36 available with HQD or HQF prefix; Size 4&6 available with HQE 8 use						
4mm	HQE4AU	2	+/- 10	N/A	+/- 6	7.2
6mm	HQE6AU	5	+/- 10	N/A	+/- 8	4.8
8mm	HQ_8AU	10	+/- 10	+/- 30	+/- 15	3.6
10mm	HQ_10AU	10	+/- 10	+/- 30	+/- 15	3.6
12mm	HQ_12AU	20	+/- 10	+/- 60	+/- 30	2.4
14mm	HQ_14AU	20	+/- 10	+/- 60	+/- 30	2.4
16mm	HQ_16AU	50	+/- 10	+/- 150	+/- 75	1.8
20mm	HQ_20AU	50	+/- 10	+/- 150	+/- 75	1.8
24mm	HQ_24AU	100	+/- 10	+/- 300	+/- 150	1.2
28mm	HQ_28AU	100	+/- 10	+/- 300	+/- 150	1.2
32mm	HQ_32AU	200	+/- 10	+/- 600	+/- 300	0.9
36mm	HQ_36AU	200	+/- 10	+/- 600	+/- 300	0.9

HQC/HQE - 8 or 16 use
 HQG - single use
 HQN - single/8/16 use
 HQD/HQF - 50 use

* HQC & HQN only available in sizes 4,6,8,10
 ** HQD also available for Japan only
 *** HQD also available for Taiwan only