



Reverse Osmosis Water Systems

Purer water for better reprocessing results

Reverse osmosis water systems are ideal for sterile processing departments; they efficiently remove contaminants including sediment, dissolved solids, and microorganisms.

Unlike water softeners that add instrument-degrading salt to your water supply, reverse osmosis systems create critical water, required for final rinse. When used throughout your process, RO water can improve your cleaning effectiveness, optimize steam quality in autoclaves, and enhance patient safety.

Easy to install and operate

- Compact space saving design
- Components easily accessible
- Pre-plumbed, wired, and assembled
- Factory tested and preserved
- Low operation costs
- Low Maintenance Costs
- Easy Maintenance and Servicing
- 1-Year Limited Warranty



CSROCT01

7000-gallon capacity system

A Reverse Osmosis Water System removes sediment, dissolved solids, microorganisms, and other impurities from your water supply for more effective instrument processing.

Learn more:

[Water for Instrument Processing](#)

Infection Control Today

Available in various sizes depending on water usage levels

All models include sediment and carbon block pre-filters, low energy membranes, and a high-pressure pump all mounted onto a corrosion-resistant, powder-coated aluminum frame.

The 7000-gallon capacity system includes A TDS monitor, multi-stage high-pressure pump, low energy membranes and three integrated pre-filters as standard features.

Product Specifications


Part ID	System Dimensions (L x W x H)	Vessels	Gals/day
CSROAT02	14 x 20 x 27 in. / 36 x 51 x 69 cm	2	500
CSROAT01	14 x 20 x 27 in. / 36 x 51 x 69 cm	3	1000
CSROBT01	19 x 23 x 46 in. / 48 x 58 x 116 cm	3	2000
CSROCT01	30 x 38 x 47 in. / 76 x 96 x 119 cm	4	7000

Available Accessories

Part ID	Description	Vessels	Gals/day
CSROTANK	Storage tank to collect RO water	2	Up to 500
CSROUV01	UV water purification system add-on	3	Up to 1000



CSROAT01
1000-gallon capacity system



Ask us for help determining the right model for your needs.

Case Medical, Inc.
50 West Street
Bloomfield, NJ 07003
U.S. and Canada 1-888-227-CASE
Worldwide +1 201-313-1999
www.casemed.com

Copyright © 2021 by Case Medical, Inc.
All rights reserved. Rev A 4/21



For more information contact
info@casemed.com
Phone: 888-227-CASE (2273)
www.casemed.com