


Your Guide to Life365's Online Portal



CARE COORDINATOR: UNDERSTANDING PORTAL BASICS


Coordinator Portal Core Functions

- 1 Review Member Populations Health Data
- 2 Establish Vital Thresholds for Members
- 3 Review Notifications for Breached Threshold Parameters
- 4 Assign, Remove, & Reassign Smart Devices
- 5 Graphing & Trending Functions
- 6 Create Reminders for Members
- 7 Review Survey Responses

Contents

Login	1
Member List Dashboard	2
Support Material Access and Account Actions.....	3
Change Portal Language	3
Change Account Password.....	3
Logout	3
Forgot Password.....	4

Login



Life365

English

Coordinator Portal

Email

Password

[Forgot your password?](#)

Remember me

A

Sign into the Coordinator Portal with a valid email & password.

Select Sign In (A).

<https://OrgName.ent.mprd.life365inc.com/en/user/login>

Member List Dashboard

The screenshot shows the Life365 Member List Dashboard. At the top left is the Life365 logo. To the right are links for 'Support', 'English', and 'Life365 Coordinator'. Below the logo are tabs for 'Vitals' and 'Notifications 12'. A search bar labeled 'Search Users' with 'Go' and 'Clear' buttons is on the right. The main content area has a blue header with 'Member' and 'Last Measurement' (with a dropdown arrow). Below this, a member profile for 'Jen Smith' is shown, including coordinator, group, site, and IMEI. To the right of the profile are several vital signs: BP (112/75), ECG (N/A), Glc (106 mg/dL), HR (N/A), SpO2 (97%), Respiration Rate (N/A), SP (N/A), Temp (97.7 °F), and Wgt (145.3 lbs).

A — **Members List**


B — **Jen Smith**
Coordinator: Emily Fellows
Group: Life365, Inc.
Site: Tempe
IMEI: 123456789012345

C — **Last Measurement**


Monitor your Member population within the Member List Dashboard **(A)**.

View your Member population by Member **(Name or Smart Device IMEI)** **(B)**, and Latest Measurement **(Most Recent Vitals)** **(C)**.

Support Material Access and Account Actions



Access Support Materials such as Bluetooth Device Pairing Quick Tip Guides, & Portal Training Materials, by selecting the Support Tab located in the top right of your screen (A).




Change your account Language by Selecting the Language option (B).

Change your account Password by Selecting the option Change Password (C).

Logout from your account by Selecting the option Logout (D).

Forgot Password



The screenshot shows the Life365 Coordinator Portal login interface. At the top left is the Life365 logo, and at the top right is a language selector set to 'English'. A blue header bar contains the text 'Coordinator Portal'. Below this, there are two text input fields: 'Email' and 'Password'. A blue link labeled 'Forgot your password?' is located below the password field, highlighted with an orange box and the letter 'A'. Below the password field is a 'Remember me' checkbox. At the bottom of the form is a green 'Sign In' button.

On the Coordinator Portal Login page, Select the [Forgot your Password? \(A\)](#) link.

Enter the email address associated with your Coordinator Portal Account.

Allow a few moments to receive the password reset email. Follow the password reset email instructions to create a new password.

If you do not receive this email, **contact your Administration** for assistance.