

# Telehealth Cabinet

Ergonomics at work to keep your equipment secure and ready.

A complete solution for your telehealth requirements, the Telehealth Cabinet keeps all your accessories along with the monitor, camera and ProEX Telehealth Hub in one secure and organised location.

All equipment is pre-wired so the cabinet need only be opened and powered-on before use. The camera and monitor are perfectly positioned to provide a user-friendly interface to the remote physician.



**Easy to install  
and prewired.**



# Telehealth Cabinet

Ergonomics at work to keep your equipment secure and ready.

Illumination is provided through an integrated LED bar in the lid of the cabinet and the fold down front door converts to the perfect work surface while the cabinet is in use. All doors are connected with gas struts to guarantee smooth and simple operation while opening and closing the cabinet.

Installation is easy with custom brackets for the wall. Once the brackets are securely mounted, the cabinet is simply hung on the brackets and locked into place.

The Telehealth Cabinet allows your clinical space to be used for multi-purpose activities safe in the knowledge that your tele-health equipment is secure and ready for use at a moments notice.

The cabinet may also be configured to facilitate unattended access to the contents by remotely unlocking the doors over the internet.

Note: Cabinet to be installed in accordance with user manual by a qualified builder/technician.



## OPTIONS

- Floor or wall mount brackets
- Key or combination lock
- Remote unlocking via Internet for unattended access scenarios
- PTZ 10x zoom camera
- Standard webcam
- 32" high definition monitor

## SPECIFICATIONS

<b>Dimensions</b>	W1350mm x H1000mm x D300mm
<b>Power</b>	110-240V AC
<b>Exposure control</b>	Automatic AGC
<b>Installation</b>	Qualified technician only
<b>Monitor</b>	Max 32" monitor
<b>Illumination</b>	Integrated LED
<b>ProEX</b>	Mounted and pre-wired

For more information or to book a demonstration, please contact Visionflex at [sales@visionflex.com](mailto:sales@visionflex.com)

+61 2 8914 4000

 [visionflex.com](http://visionflex.com)



**visionflex**  
Equality healthcare