

The POINT DIGIT II

by point designs



point designs

www.pointdesignsllc.com | sales@pointdesignsllc.com

POINT DIGIT II

MCP center
of rotation

Titanium construction,
150 lb load capacity

Gunmetal or
polished finish

Release button &
auto spring back

Touchscreen compatible
fingertip pads for enhanced grip

Ratcheting
positioning
mechanism

Anatomical
flexion

MOUNTING KIT

Torx® mounting screws

Steel construction

Breakaway design for
1-4 finger fittings

Includes lamination
spacers

6 stock sizes with fast delivery



point designs

Thank you for your interest in the Point Digit II. Our team at Point Designs is committed to the innovation of devices and algorithms in the upper limb prosthetic field. We develop products that provide prosthetists and users with solutions to meet the unique needs of every patient. The Point Digit II provides an effective and robust solution for many upper limb amputees with finger loss and provides your clients with the capability to enjoy their work, family, and life to the fullest.