



## Beauty and Precision



**PATIENTS** today are better informed and demand beauty and precision. The expectation of patients is increased functionality with esthetically pleasing restorations at a competitive price. Patients want great, boutique laboratory restorations, but are often not willing to accept the cost associated with these restorations.

**PROLINE** is the answer for the demanding and price conscious patient. Boutique quality at competitive prices from a laboratory with over 40 years of restorative experience, with the quality and service you expect. Reliable, beautiful restorations that increase patient acceptance and keep costs manageable.

**ADVANTAGES** - Each PROLINE case is assigned to our PROLINE team, comprised of technicians certified and trained in industry leading programs such as Dawson, LVI, and CDT. Each case will be handled by this dedicated team, and will undergo eight different quality control checkpoints during production.



## PROLINE REMOVABLES

Each of our PROLINE dentures are hand crafted by expert technicians with years of experience creating beautiful removable prosthetics.

Features:

- Extended 2 year warranty
- Ivoclar Processing
- Ivoclar Phonares II Teeth
- PROLINE Custom Setup
- Custom Gingival Tinting



## PROLINE CROWN AND BRIDGE

A crown is more than just a tooth repair. It's a restoration that gives your patient confidence and a natural looking smile. For patients who want advanced characterization or a highly esthetic restoration, we provide the answer with PROLINE crown and bridge.

Features:

- Extended 7 year warranty
- All premium materials available
- Advanced Shading
- Advanced Characterization

2355 Centennial Circle, Gainesville, GA 30504

[www.pittmandental.com](http://www.pittmandental.com) | 800-235-4720





# PROLINE Restorations Fee Schedule

## Fixed and Removables Restorations

### All Ceramic Crown and Bridge

#### Diamante™ Zirconia

Crown.....	231
Bridge.....	254

Our Diamante™ Zirconia bruxer is ideal for posterior restorations; especially bruxer patients. We mill from a medical grade zirconia and use a proprietary shading process for superior esthetics. Indications: Posterior, Bruxism patients.

#### Diamante™ Layered Zirconia

Crown.....	240
Bridge.....	266

We take the same imported premium grade zirconia we use for our Diamante™ bruxers and stack multiple layers of porcelain to create an esthetic and translucent crown capable of rivaling all leading brands. Indications: Anterior, Posterior, Esthetic Zone.

#### LuMax™

Crown.....	231
Crown Layered.....	254
Inlay/Onlay.....	231

LuMax is a high translucent zirconia restoration that offers the esthetics of traditional e.max® with increased flexural strength. Indications: Anterior, Posterior, Veneers, Inlay/Onlay, Esthetic Zone.

#### IPS e.max®

Crown/Veneer Pressed.....	231
Crown Layered.....	254
Thin Press.....	254

IPS e.max® uses a patented lithium disilicate formulation optimizing CAD/CAM and pressing design processes with a high performance ceramic material. Indications: Anterior, Posterior, Veneers, Inlay/Onlay, Esthetic Zone.

### Porcelain Fused to Metal

Non-Precious.....	184
Palladium Nobel Alloy <sup>1</sup> .....	201
56% White Gold <sup>2</sup> .....	211
76% Yellow Gold <sup>2</sup> .....	222

Porcelain to Gold restorations include 1 D.W.T. of alloy.

<sup>1</sup>Palladium Alloys are based on market value of \$400, for each \$50 increase over \$400, add \$2 per unit.

<sup>2</sup>Gold prices are based on a market value of \$400 an ounce. For each \$50 increase over the base of \$400, add \$5 per unit.

### Dentures

Complete PROLINE Denture.....	305
-------------------------------	-----

Includes: Models, Articulation, Ivoclar Phoneras II Teeth, Custom Gingival Tinting and Ivo Base Processing

### PROLINE Smile Kit

This kit allows the clinician to improve accuracy and efficiency of full dentures, while ensuring patient satisfaction. Each component of the kit is designed to address specific steps in diagnosis and prescription procedure. Contact us today for more information.

## Recommended Steps for a PROLINE Restoration

To ensure a successful result, our PROLINE team recommends submitting the following with each PROLINE case:

Complete Treatment Plan, Patient Expectations and Esthetic Requests, Provisional Model or Study Model, Full Face Smile and Retracted Smile Digital Pictures, Full Arch/Stock Tray Impressions or Digital Impression, Separate Bite Registration(s) Over Preps, and Smile Design Selection (if applicable). Please feel free to contact us with any questions.

Support@pittmandental.com

2355 Centennial Circle, Gainesville, GA 30504

www.pittmandental.com | 800-235-4720

