

RAW
DATA

RAW DATA MADE EASY

MUNICH - BERLIN - ROSTOCK - PHOENIX - WASHINGTON, D.C.



BM-Viewer



BM-Viewer is an agent that is able to monitor the file storage of a laboratory device. New or modified files can be transferred to the BM-Viewer ECM. When changes are made, new versions of the file are automatically created.

BM-Viewer is also capable of

automatically retrieving metadata, e.g. from file paths, text files or files in proprietary formats and databases. Reasons for changes can also be requested from the user. In BM-Viewer all files are provided with an audit trail, changes can be traced via the versioning.

- ✓ Access control
- ✓ Audit trails
- ✓ File versions
- ✓ Electronic signature
- ✓ Lifecycle
- ✓ Automated indexing
- ✓ E- Long term archive
- ✓ Instrument monitoring

THE MONITORING STRATEGY

BM-Viewer can be configured in different ways. One possibility is the immediate transfer of the data to the archive, but alternatively the transfer can also take place after a certain time or after the creation of a certain file. In the latest version, the BM-Viewer can now also be configured and monitored via a central




server. The configuration is done with the help of modules that provide functionality such as the determination of metadata. These modules are configured using XML files. All configuration changes are traceable in a system audit trail. The configuration server supports device monitoring via a status page. The current status is visible at a glance; systems in error state are distinguished from correctly closed systems.



Based on windream

BM-windream is the 21 CFR Part 11 compliant data and document management system for the pharmaceutical industry. It is based on the market leading windream document management system from windream GmbH, Bochum. In addition to our industry-specific modules, we have access to a series of tried and tested windream add-on modules from our partner network, e.g. for e-mail archiving, linking to ERP systems, CRM systems, knowledge management or scanning solutions. BM-windream provides a basis for your data and document management that can be validated at relatively low cost (configurable standard software complying with GAMP5) and is easy to use. Users access your documents using a virtual network drive in Windows Explorer, which means that little training is required. BM-windream provides a standardised user interface over the entire life cycle of a file, from creation to deletion. At all phases in this life cycle, a change history (audit trail) is recorded.



CONTACT

EU
P: +49 30 77018110
E: info@arivis.com

BioMedion GmbH
Am Borsigturm 100,
13507 Berlin, Germany

US
P: +1 602 957 2160
E: info@arivis.com

arivis, Inc.
3514 N. Power Road, Suite 115
Mesa, AZ 85215, USA

MUNICH - BERLIN - ROSTOCK - PHOENIX - WASHINGTON, D.C.